



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116472

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3868-05	33458	1	1.00	2.00	2.00	100.0	OIL SAMPLE
M3892-01	S-829-J-SO-1-1.5-092221	2	1.15	9.78	7.84	77.5	
M3892-02	S-829-J-SO-1-1.5-092221 -FD	3	1.14	9.96	8.2	80.0	
M3892-03	S-829-J-SO-2-2.5-092221	4	1.16	9.91	7.9	77.0	
M3892-04	S-829-J-SO-3-3.5-092221	5	1.15	9.94	8.37	82.1	
M3892-05	S-822-N-SO-0-0.5-092221	6	1.16	9.68	6.53	63.0	
M3892-06	S-822-N-SO-1-1.5-092221	7	1.12	9.87	8.54	84.8	
M3892-07	S-823-N-SO-0-0.5-092221	8	1.19	9.81	8.89	89.3	
M3892-08	S-823-N-SO-1-1.5-092221	9	1.16	9.96	9.06	89.8	
M3892-09	S-824-N-SO-0-0.5-092221	10	1.16	9.96	8.11	79.0	
M3892-10	S-824-N-SO-1-1.5-092221	11	1.18	9.72	8.13	81.4	
M3892-11	S-825-N-SO-0-0.5-092221	12	1.16	9.93	8.31	81.5	
M3892-12	S-825-N-SO-1-1.5-092221	13	1.14	9.76	8.29	82.9	
M3892-13	S-826-N-SO-0-0.5-092221	14	1.13	9.71	9.25	94.6	
M3892-14	S-826-N-SO-1-1.5-092221	15	1.15	9.98	8.1	78.7	
M3892-15	S-832-J-SO-1-1.5-092221	16	1.12	9.88	8.4	83.1	
M3892-16	M3892-15MS	17	1.12	9.88	8.4	83.1	
M3892-17	M3892-15MSD	18	1.12	9.88	8.4	83.1	
M3892-18	S-832-J-SO-2-2.5-092221	19	1.18	9.95	7.84	75.9	
M3892-19	S-832-J-SO-2-2.5-092221 -FD	20	1.18	9.97	8.2	79.9	
M3892-20	S-832-J-SO-3-3.5-092221	21	1.13	9.95	8.63	85.0	
M3899-01	PCB-01	23	1.14	9.56	8.7	89.8	
M3899-02	PCB-02	24	1.19	9.65	8.65	88.2	
M3899-03	DIESEL-SOIL-01	25	1.16	9.98	8.99	88.8	
M3903-01	225	26	1.14	9.89	9.11	91.1	
M3903-03	RBR-254	27	1.15	9.95	9.13	90.7	
M3903-05	11976	28	1.13	9.71	9.36	95.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116472

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3903-07	202	29	1.16	9.95	8.9	88.1	
M3904-01	JET-GROUT-SOIL-STEEL-WH EEL	22	1.13	9.72	6.74	65.3	
M3906-01	BS-DRUM-COMP	30	1.00	2.00	2.00	100.0	oily-stone
M3910-01	11986	31	1.14	9.94	9.56	95.7	
M3910-03	208	32	1.17	9.56	8.9	92.1	
M3914-01	LV-DRUM-SOIL	33	1.14	9.63	9.23	95.3	
M3914-03	LV-DRUM-PADS	34	1.00	2.00	2.00	100.0	DRUM-PADS
M3916-01	CONC-1	35	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-02	CONC-1-EPH-2	36	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-04	CONC-2	37	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-05	CONC-2-EPH-2	38	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-07	CONC-3	39	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-08	CONC-3-EPH-2	40	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-10	CONC-4	41	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-11	CONC-4-EPH-2	42	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-13	BRICK-1	43	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3916-14	BRICK-1-EPH-2	44	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3917-01	TP-3	45	1.16	9.81	8.69	87.1	
M3917-02	TP-3-EPH-2	46	1.16	9.81	8.69	87.1	
M3917-04	TP-6	47	1.11	9.7	8.48	85.8	
M3917-05	TP-6-EPH-2	48	1.11	9.7	8.48	85.8	
M3917-07	TP-10	49	1.14	9.94	8.97	89.0	
M3917-08	TP-10-EPH-2	50	1.14	9.94	8.97	89.0	
M3917-10	TP-9	51	1.16	9.72	8.66	87.6	
M3917-11	TP-9-EPH-2	52	1.16	9.72	8.66	87.6	
M3917-13	TP-8	53	1.13	9.84	8.88	89.0	
M3917-14	TP-8-EPH-2	54	1.13	9.84	8.88	89.0	
M3917-16	TP-4	55	1.19	9.88	8.38	82.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116472

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3917-17	TP-4-EPH-2	56	1.19	9.88	8.38	82.7	
M3917-19	TP-1	57	1.1	9.76	8.44	84.8	
M3917-20	TP-1-EPH-2	58	1.1	9.76	8.44	84.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$