



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 9/25/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:30
In Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:55
Out Date: 09/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116490

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3852-01	EX367	1	1.14	9.71	8.75	88.8	
M3852-02	EX368	2	1.12	9.8	8.91	89.7	
M3852-03	EX369	3	1.16	9.76	8.79	88.7	
M3852-04	EX370	4	1.13	9.71	8.4	84.7	
M3852-05	EX370D	5	1.13	9.71	8.4	84.7	
M3852-06	EX370S	6	1.13	9.71	8.4	84.7	
M3852-07	EX371	7	1.15	9.54	8.3	85.2	
M3852-08	EX372	8	1.17	9.68	8.26	83.3	
M3852-09	EX373	9	1.11	9.86	8.58	85.4	
M3852-10	EX374	10	1.15	9.5	8.3	85.6	
M3852-11	EX375	11	1.16	9.97	8.89	87.7	
M3852-12	EX376	12	1.14	9.7	7.88	78.7	
M3852-13	EX377	13	1.16	9.81	8.26	82.1	
M3852-14	EX378	14	1.15	9.96	8.4	82.3	
M3852-15	EX379	15	1.15	9.71	8.39	84.6	
M3852-16	EX380	16	1.16	9.78	8.5	85.2	
M3852-17	EX381	17	1.19	9.74	8.15	81.4	
M3852-18	EX382	18	1.14	9.98	8.94	88.2	
M3852-19	EX383	19	1.1	9.65	8.5	86.5	
M3852-20	EX385	20	1.16	9.55	8.1	82.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$