



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/27/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 11:00
In Date: 09/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-01

OVENTEMP OUT Celsius(°C): 103
Time OUT: 11:10
Out Date: 09/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116504

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3905-16	X3725	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-17	X3725D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-18	X3725S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-19	X3726	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-20	X3726D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-21	X3726S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-22	X3727	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-23	X3727D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-24	X3727S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-25	X3728	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-26	X3728D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3905-27	X3728S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$