



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/27/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:15
In Date: 09/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116510

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3877-01	EW2L6	1	1.15	9.71	8.89	90.4	
M3877-02	EW2L7	2	1.13	9.51	8.68	90.1	
M3877-03	EW2L8	3	1.16	9.78	8.57	86.0	
M3877-04	EW2L9	4	1.19	9.7	8.95	91.2	
M3877-05	EW2M1	5	1.16	9.88	8.89	88.6	
M3877-06	EW2M2	6	1.19	9.61	8.17	82.9	
M3877-07	EW2M3	7	1.12	9.94	8.98	89.1	
M3877-08	EW2M4	8	1.19	9.78	9.24	93.7	
M3877-09	EW2M5	9	1.18	9.95	9.3	92.6	
M3877-10	EW2M6	10	1.13	9.63	8.82	90.5	
M3877-11	EW2M7	11	1.16	9.64	9.11	93.7	
M3877-12	EW2M8	12	1.16	9.88	9.05	90.5	
M3877-13	EW2M9	13	1.17	9.57	8.94	92.5	
M3877-14	EW2N0	14	1.17	9.6	8.75	89.9	
M3877-15	VHBLK001	15	1.00	2.00	2.00	100.0	VHBLK
M3877-16	EW2L8MS	16	1.16	9.78	8.57	86.0	
M3877-17	EW2L8MSD	17	1.16	9.78	8.57	86.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$