



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 10/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116852

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4242-09	T1-4-BOTTOM-GRAB	1	1.15	9.8	9.39	95.3	
M4242-10	T1-COMP- (1-4)	2	1.16	9.75	9.28	94.5	
M4242-20	T2-4-BOTTOM-GRAB	3	1.17	9.67	9.03	92.5	
M4242-21	T2-COMP- (1-4)	4	1.12	9.54	9.08	94.5	
M4243-09	T3-4-BOTTOM-GRAB	5	1.16	9.85	9.00	90.2	
M4243-10	T3-COMP- (1-4)	6	1.15	9.79	8.66	86.9	
M4243-20	T4-4-BOTTOM-GRAB	7	1.16	9.76	8.71	87.8	
M4243-21	T4-COMP- (1-4)	8	1.17	9.95	9.48	94.6	
M4244-09	T5-4-BOTTOM-GRAB	9	1.16	9.94	9.00	89.3	
M4244-10	T5-COMP- (1-4)	10	1.16	9.62	8.23	83.6	
M4248-01	CHVB-0513	29	1.15	9.98	1.24	1.0	GEL MATRIX
M4248-02	CHVB-0513-EPH-2	30	1.15	9.98	1.24	1.0	GEL MATRIX
M4248-03	72-12002	31	1.16	9.94	9.59	96.0	
M4248-04	72-12002-EPH-2	32	1.16	9.94	9.59	96.0	
M4249-01	TP1- (S)	11	1.12	9.77	8.83	89.1	
M4249-03	TP1- (S)	12	1.12	9.77	8.83	89.1	
M4249-04	TP1- (D)	13	1.16	9.63	8.74	89.5	
M4249-06	TP1- (D)	14	1.16	9.63	8.74	89.5	
M4249-07	TP2- (S)	15	1.19	9.83	9.17	92.4	
M4249-09	TP2- (S)	16	1.19	9.83	9.17	92.4	
M4249-10	TP2- (D)	17	1.11	9.62	8.67	88.8	
M4249-12	TP2- (D)	18	1.11	9.62	8.67	88.8	
M4250-01	PIPESEG-A	19	1.00	2.00	2.00	100.0	WIPE SAMPLE
M4250-02	PIPESEG-B	20	1.00	2.00	2.00	100.0	WIPE SAMPLE
M4252-01	TEST-BA- (1)	21	1.12	9.82	9.06	91.3	
M4252-03	TEST-BA- (1)	22	1.12	9.82	9.06	91.3	
M4252-04	TEST-BB- (1)	23	1.15	9.95	9.38	93.5	
M4252-06	TEST-BB- (1)	24	1.15	9.95	9.38	93.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 10/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116852

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4252-07	MH-CA- (1)	25	1.16	9.98	9.42	93.7	
M4252-09	MH-CA- (1)	26	1.16	9.98	9.42	93.7	
M4252-10	MH-CB- (1)	27	1.11	9.73	8.42	84.8	
M4252-12	MH-CB- (1)	28	1.11	9.73	8.42	84.8	
M4253-01	131	33	1.00	2.00	2.00	100.0	oily rags
M4254-05	SVOC-GPC-BLANK	34	1.00	2.00	2.00	100.0	
M4254-06	PEST-GPC-BLANK	35	1.00	2.00	2.00	100.0	
M4254-07	PEST-GPC-BLANK-SPIKE	36	1.00	2.00	2.00	100.0	
M4254-08	PCB-GPC-BLANK	37	1.00	2.00	2.00	100.0	
M4254-09	PCB-GPC-BLANK-SPIKE	38	1.00	2.00	2.00	100.0	
M4254-10	SVOC-GPC2-BLANK	39	1.00	2.00	2.00	100.0	
M4254-11	PEST-GPC2-BLANK	40	1.00	2.00	2.00	100.0	
M4254-12	PEST-GPC2-BLANK-SPIKE	41	1.00	2.00	2.00	100.0	
M4254-13	PCB-GCP2-BLANK	42	1.00	2.00	2.00	100.0	
M4254-14	PCB-GCP2-BLANK-SPIKE	43	1.00	2.00	2.00	100.0	
M4256-01	454BASIN-20211014	44	1.11	9.88	9.55	96.2	
M4256-02	454BASIN-20211014	45	1.11	9.98	9.55	95.2	
M4259-01	WC-1-CS	46	1.13	9.52	8.78	91.2	
M4259-03	WC-1-CS-EPH-2	47	1.13	9.52	8.78	91.2	
M4259-04	WC-1-BS	48	1.18	9.95	9.4	93.7	
M4259-06	WC-1-BS-EPH-2	49	1.18	9.95	9.4	93.7	
M4259-07	WC-1-BD	50	1.13	9.71	9.07	92.5	
M4259-09	WC-1-BD-EPH-2	51	1.13	9.71	9.07	92.5	
M4259-10	WC-1-AS	52	1.14	9.56	8.83	91.3	
M4259-12	WC-1-AS-EPH-2	53	1.14	9.56	8.83	91.3	
M4260-01	TEST-DA- (1)	54	1.16	9.84	8.68	86.6	
M4260-03	TEST-DA- (1)	55	1.16	9.84	8.68	86.6	
M4260-04	TEST-DB- (1)	56	1.17	9.97	9.12	90.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 10/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 10/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 10/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116852

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4260-06	TEST-DB- (1)	57	1.17	9.97	9.12	90.3	
M4261-01	S-819-N-S-1-1.5-093021	58	1.15	9.91	9.24	92.4	
M4261-02	S-819-N-S-1-1.5-093021-FD	59	1.17	9.63	8.89	91.3	
M4261-03	S-818-N-SO-1-1.5-100121	60	1.14	9.82	9.37	94.8	
M4261-04	S-818-N-SO-3-3.5-100121	61	1.16	9.86	8.39	83.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$