

Analysis Method: 9034
 Parameter: Reactive Sulfide
 Run Number: LB116949

ANALYST: ketankumar
 SUPERVISOR REVIEW BY: apatel
 Constant: 16000
 Normality1: 0.025
 Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W2843
IODINE SOLUTION .025N 1L	W2864
Starch Solution, 4L	W2703

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB140080BL		1	5.00	50	2.00	0.00	1.90	1.90	0.10	0.00	0.00	10/21/2021	09:15
2	M4269-03		1	5.03	50	2.00	1.90	3.74	1.84	0.16	0.06	4.77	10/21/2021	09:18
3	M4269-03DUP		1	5.04	50	2.00	3.74	5.60	1.86	0.14	0.04	3.17	10/21/2021	09:20
4	M4269-05		1	5.05	50	2.00	5.60	7.42	1.82	0.18	0.08	6.34	10/21/2021	09:22
5	M4269-07		1	5.07	50	2.00	7.42	9.30	1.88	0.12	0.02	1.58	10/21/2021	09:25
6	M4269-09		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.04	3.16	10/21/2021	09:28
7	M4269-11		1	5.03	50	2.00	1.86	3.70	1.84	0.16	0.06	4.77	10/21/2021	09:30
8	M4269-13		1	5.07	50	2.00	3.70	5.52	1.82	0.18	0.08	6.31	10/21/2021	09:32
9	M4269-15		1	5.09	50	2.00	5.52	7.38	1.86	0.14	0.04	3.14	10/21/2021	09:35
10	M4274-02		1	5.08	50	2.00	0.00	1.88	1.88	0.12	0.02	1.57	10/21/2021	09:38
11	M4274-05		1	5.07	50	2.00	1.88	3.72	1.84	0.16	0.06	4.73	10/21/2021	09:40
12	M4274-08		1	5.06	50	2.00	3.72	5.60	1.88	0.12	0.02	1.58	10/21/2021	09:42

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume