



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 11/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117180

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4477-01	T1-EP01-N	1	1.14	8.67	9.81	9.36	94.8	
M4477-02	T1-EP02-NE	2	1.12	8.56	9.68	9.33	95.9	
M4477-03	T1-EP03-S	3	1.15	8.49	9.64	9.21	94.9	
M4477-04	T1-B1-S-N	4	1.1	8.51	9.61	9.1	94.0	
M4477-05	T1-B2-NW	5	1.13	8.73	9.86	9.34	94.0	
M4477-06	T5-EP1-S	6	1.15	8.81	9.96	9.42	93.9	
M4477-07	T5-EP2-NW	7	1.13	8.44	9.57	8.76	90.4	
M4477-08	T1-B3-N	8	1.13	8.73	9.86	9.44	95.2	
M4477-09	T5-EP3-N	9	1.16	8.44	9.6	8.98	92.7	
M4477-10	T5-B1-NE	10	1.13	8.55	9.68	9.06	92.7	
M4477-11	T5-B2-NW	11	1.11	8.77	9.88	9.37	94.2	
M4485-01	AG-2	12	1.14	8.81	9.95	9.47	94.6	
M4485-02	AG-3	13	1.13	8.73	9.86	9.23	92.8	
M4485-03	AG-4	14	1.12	8.62	9.74	9.26	94.4	
M4485-04	AG-5	15	1.15	8.58	9.73	8.98	91.3	
M4485-05	AG-6	16	1.18	8.79	9.97	9.16	90.8	
M4485-06	AG-7	17	1.13	8.68	9.81	9.13	92.2	
M4485-07	AG-8	18	1.18	8.42	9.6	9.03	93.2	
M4485-08	AG-9	19	1.18	8.46	9.64	8.66	88.4	
M4485-09	AG-10	20	1.18	8.45	9.63	8.87	91.0	
M4497-07	346	21	1.14	8.39	9.53	8.6	88.9	
M4519-01	SB-1	22	1.14	8.56	9.7	9.06	92.5	
M4519-02	SB-2	23	1.16	8.79	9.95	8.44	82.8	
M4519-03	SB-3	24	1.16	8.34	9.5	8.83	92.0	
M4519-04	SB-4	25	1.16	8.77	9.93	9.04	89.9	
M4519-05	SB-5	26	1.15	8.69	9.84	9.06	91.0	
M4520-01	40540	27	1.15	8.60	9.75	8.93	90.5	
M4523-01	VWE-B4-110421-0-10	28	1.15	8.36	9.51	8.82	91.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 11/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117180

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4523-02	VWE-B4-110421-10-20	29	1.15	8.15	9.3	8.32	88.0	
M4523-03	VWE-B4-110421-20-30	30	1.15	8.54	9.69	8.45	85.5	
M4523-04	VWE-B4-110421-0-1	31	1.15	8.62	9.77	9.08	92.0	
M4523-05	VWE-B4-110421-14-15	32	1.15	8.70	9.85	8.89	89.0	
M4523-06	VWE-B4-110421-29-30	33	1.15	8.45	9.6	8.87	91.4	
M4523-09	VWE-B4-110421-BD	34	1.15	8.65	9.8	8.43	84.2	
M4523-10	VWE-B4-110421-BD	35	1.15	8.82	9.97	9.24	91.7	
M4525-01	SEC-2021-26A	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-02	SEC-2021-26B	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-03	SEC-2021-27A	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-04	SEC-2021-27B	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-05	SEC-2021-28A	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-06	SEC-2021-28B	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-07	SEC-2021-29A	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-08	SEC-2021-29B	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-09	SEC-2021-30A	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-10	SEC-2021-30B	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-11	SEC-2021-31A	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-12	SEC-2021-31B	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-13	SEC-2021-32A	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4525-14	SEC-2021-32B	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4529-01	6-INCH-STEEL	50	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
M4529-02	12-INCH-STEEL	51	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
M4536-01	SEC-2021-33A	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-02	SEC-2021-33B	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-03	SEC-2021-34A	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-04	SEC-2021-34B	55	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-05	SEC-2021-35A	56	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 11/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117180

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4536-06	SEC-2021-35B	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-07	SEC-2021-36B	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-08	SEC-2021-37A	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-09	SEC-2021-38B	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-10	SEC-2021-39A	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4536-11	SEC-2021-39B	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4537-01	S1-S2-S3	81	1.12	8.41	9.53	8.16	83.7	
M4538-01	WC-C1-COMP	63	1.15	8.60	9.75	8.76	88.5	
M4539-01	WC-C2-COMP	64	1.15	8.66	9.81	8.84	88.8	
M4539-03	M4539-01MS	65	1.15	8.66	9.81	8.84	88.8	
M4539-04	M4539-01MSD	66	1.15	8.66	9.81	8.84	88.8	
M4540-01	TP-1	67	1.15	8.39	9.54	8.48	87.4	
M4540-02	TP-1-EPH-2	68	1.15	8.45	9.6	8.84	91.0	
M4540-03	TP-2	69	1.15	8.46	9.61	8.29	84.4	
M4540-04	TP-2-EPH-2	70	1.15	8.61	9.76	8.4	84.2	
M4540-05	TP-3	71	1.15	8.70	9.85	8.59	85.5	
M4540-06	TP-3-EPH-2	72	1.15	8.37	9.52	7.8	79.5	
M4540-07	TP-4	73	1.15	8.45	9.6	8.65	88.8	
M4540-08	TP-4-EPH-2	74	1.15	8.59	9.74	8.79	88.9	
M4540-09	TP-5	75	1.15	8.45	9.6	8.67	89.0	
M4540-10	TP-5-EPH-2	76	1.15	8.79	9.94	8.91	88.3	
M4540-11	TP-6	77	1.15	8.70	9.85	8.09	79.8	
M4540-12	TP-6-EPH-2	78	1.15	8.72	9.87	8.19	80.7	
M4541-01	WC-3-AS	79	1.15	8.73	9.88	9.06	90.6	
M4541-03	WC-3-AS-EPH-2	80	1.15	8.72	9.87	8.74	87.0	
M4544-01	MOO-2021-10A	98	1.00	1.00	2.00	2.00	100.0	PILC
M4544-02	MOO-2021-10B	99	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4544-03	MOO-2021-10C	100	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 11/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117180

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4544-04	MOO-2021-11A	101	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4544-05	MOO-2021-11B	102	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4544-06	MOO-2021-11C	103	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4544-07	MOO-2021-11D	104	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4548-01	LAW-2021-1A	105	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4548-02	LAW-2021-2A	106	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4548-03	LAW-2021-2B	107	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4548-04	LAW-2021-3A	108	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4548-05	LAW-2021-3B	109	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4561-01	ETGI-280	82	1.00	1.00	2.00	2.00	100.0	WOOD SAMPLE
M4562-01	T-1	83	1.18	8.50	9.68	7.92	79.3	
M4562-02	T-1A	84	1.16	8.56	9.72	8.18	82.0	
M4563-01	S-1	85	1.13	8.55	9.68	8.37	84.7	
M4569-01	CL-02-110521	86	1.14	8.80	9.94	9.23	91.9	
M4569-02	CL-02-110521-EPH-2	87	1.15	8.80	9.95	9.13	90.7	
M4570-01	SOM-2021-1A	110	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-02	SOM-2021-1B	111	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-03	SOM-2021-2A	112	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-04	SOM-2021-2B	113	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-05	SOM-2021-3A	114	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-06	SOM-2021-3B	115	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-07	SOM-2021-4A	116	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-08	SOM-2021-5A	117	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-09	SOM-2021-5B	118	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-10	SOM-2021-6A	119	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-11	SOM-2021-6B	120	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-12	KMA79SC-A	121	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4570-13	KMA79SC-B	122	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/8/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 11/05/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/06/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117180

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4574-01	VWE-B6-110521-0-10	88	1.12	8.76	9.88	9.06	90.6	
M4574-02	VWE-B6-110521-10-20	89	1.1	8.69	9.79	9.22	93.4	
M4574-03	VWE-B6-110521-20-30	90	1.17	8.73	9.9	9.45	94.8	
M4574-04	VWE-B6-110521-0-1	91	1.14	8.59	9.73	9.21	93.9	
M4574-05	VWE-B6-110521-14-15	92	1.13	8.50	9.63	8.75	89.6	
M4574-06	VWE-B6-110521-29-30	93	1.13	8.61	9.74	9.23	94.1	
M4576-01	SB-1- (12-12.5)	94	1.13	8.63	9.76	8.98	91.0	
M4576-03	SB-2- (9-9.5)	95	1.16	8.44	9.6	9.44	98.1	
M4576-04	DW-2	96	1.16	8.81	9.97	7.78	75.1	
M4576-05	DW-3	97	1.18	8.60	9.78	5.81	53.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$