



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 11/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:29
Out Date: 11/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117241

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4468-01	EXAW2	1	1.13	8.44	9.57	7.57	76.3	
M4468-02	EXAW3	2	1.16	8.59	9.75	7.92	78.7	
M4468-03	EXAW4	3	1.1	8.61	9.71	7.93	79.3	
M4468-04	EXAW5	4	1.18	8.41	9.59	7.93	80.3	
M4468-05	EXAW6	5	1.18	8.40	9.58	7.39	73.9	
M4468-06	EXAW7	6	1.14	8.81	9.95	8.00	77.9	
M4468-07	EXAW8	7	1.19	8.66	9.85	8.13	80.1	
M4468-08	EXAW9	8	1.11	8.82	9.93	8.11	79.4	
M4468-09	EXAX0	9	1.13	8.37	9.5	7.25	73.1	
M4468-10	EXAX1	10	1.16	8.67	9.83	7.57	73.9	
M4468-11	EXAX2	11	1.16	8.81	9.97	8.52	83.5	
M4468-12	EXAX3	12	1.18	8.62	9.8	7.93	78.3	
M4468-13	EXAX4	13	1.12	8.69	9.81	7.99	79.1	
M4468-14	EXAX5	14	1.16	8.36	9.52	7.34	73.9	
M4468-15	EXAX6	15	1.16	8.59	9.75	7.66	75.7	
M4468-16	EXAX7	16	1.18	8.79	9.97	8.44	82.6	
M4468-17	EXAX8	17	1.19	8.55	9.74	8.23	82.3	
M4468-18	EXAX9	18	1.18	8.76	9.94	7.74	74.9	
M4468-19	EXAY0	19	1.17	8.77	9.94	8.21	80.3	
M4468-20	EXAY1	20	1.19	8.66	9.85	7.93	77.8	
M4468-21	EXAY1D	21	1.19	8.66	9.85	7.93	77.8	
M4468-22	EXAY1S	22	1.19	8.66	9.85	7.93	77.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$