



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:15
In Date: 11/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 11/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117256

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4483-01	GB815	1	1.12	8.38	9.5	9.22	96.7	
M4483-02	GB816	2	1.19	8.67	9.86	9.74	98.6	
M4483-03	GB817	3	1.16	8.50	9.66	9.51	98.2	
M4483-04	GB818	4	1.16	8.48	9.64	9.37	96.8	
M4483-05	GB819	5	1.14	8.39	9.53	9.38	98.2	
M4483-06	GB820	6	1.15	8.82	9.97	9.77	97.7	
M4483-07	GB820D	7	1.15	8.82	9.97	9.77	97.7	
M4483-08	GB820S	8	1.15	8.82	9.97	9.77	97.7	
M4483-09	GB821	9	1.16	8.54	9.7	9.42	96.7	
M4483-10	GB822	10	1.16	8.78	9.94	9.42	94.1	
M4483-11	GB823	11	1.1	8.77	9.87	9.69	97.9	
M4483-12	GB824	12	1.18	8.70	9.88	9.68	97.7	
M4483-13	GB825	13	1.16	8.82	9.98	9.39	93.3	
M4483-14	GB826	14	1.17	8.34	9.51	8.93	93.0	
M4483-15	GB827	15	1.15	8.53	9.68	9.67	99.9	
M4483-16	GB828	16	1.17	8.43	9.6	9.28	96.2	
M4483-17	GB829	17	1.17	8.80	9.97	9.71	97.0	
M4483-18	GB830	18	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$