



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:30
In Date: 11/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117257

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4486-01	GB7X3	1	1.14	8.37	9.51	6.76	67.1	
M4486-02	GB7X4	2	1.13	8.75	9.88	7.3	70.5	
M4486-03	GB7X5	3	1.15	8.63	9.78	8.24	82.2	
M4486-04	GB7X6	4	1.18	8.50	9.68	8.16	82.1	
M4486-05	GB7X7	5	1.12	8.54	9.66	7.08	69.8	
M4486-06	GB7X8	6	1.19	8.59	9.78	7.2	70.0	
M4486-07	GB7X9	7	1.15	8.82	9.97	7.81	75.5	
M4486-08	GB7Y0	8	1.14	8.80	9.94	9.17	91.3	
M4486-09	GB7Y1	9	1.18	8.79	9.97	7.52	72.1	
M4486-10	GB7Y2	10	1.15	8.81	9.96	8.36	81.8	
M4486-11	GB7Y3	11	1.16	8.71	9.87	6.23	58.2	
M4486-12	GB7Y4	12	1.13	8.39	9.52	5.98	57.8	
M4486-13	GB7Y5	13	1.14	8.63	9.77	6.55	62.7	
M4486-14	GB7Y6	14	1.18	8.80	9.98	6.47	60.1	
M4486-15	GB7Y7	15	1.15	8.81	9.96	6.69	62.9	
M4486-16	GB7Y7D	16	1.15	8.81	9.96	6.69	62.9	
M4486-17	GB7Y7S	17	1.15	8.81	9.96	6.69	62.9	
M4486-18	GB7Y8	18	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE
M4486-19	GB7Y9	19	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$