



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:42
In Date: 11/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 11/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117261

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4488-01	GB7Z0	1	1.14	8.36	9.5	9.19	96.3	
M4488-02	GB7Z1	2	1.16	8.70	9.86	9.66	97.7	
M4488-03	GB7Z2	3	1.16	8.61	9.77	9.58	97.8	
M4488-04	GB7Z3	4	1.13	8.42	9.55	9.47	99.0	
M4488-05	GB7Z4	5	1.12	8.84	9.96	9.83	98.5	
M4488-06	GB7Z5	6	1.13	8.37	9.5	9.26	97.1	
M4488-07	GB7Z6	7	1.13	8.57	9.7	9.55	98.2	
M4488-08	GB7Z7	8	1.11	8.72	9.83	9.66	98.1	
M4488-09	GB7Z8	9	1.13	8.82	9.95	9.66	96.7	
M4488-10	GB7Z9	10	1.13	8.54	9.67	9.62	99.4	
M4488-11	GB800	11	1.19	8.57	9.76	9.7	99.3	
M4488-12	GB801	12	1.16	8.56	9.72	9.52	97.7	
M4488-13	GB802	13	1.14	8.76	9.9	9.58	96.3	
M4488-14	GB803	14	1.17	8.76	9.93	9.68	97.1	
M4488-15	GB803D	15	1.17	8.76	9.93	9.68	97.1	
M4488-16	GB803S	16	1.17	8.76	9.93	9.68	97.1	
M4488-17	GB804	17	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$