



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:10
In Date: 11/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 11/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117290

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4560-01	GB8C7	1	1.18	8.63	9.81	9.63	97.9	
M4560-02	GB8C8	2	1.16	8.69	9.85	9.66	97.8	
M4560-03	GB8C9	3	1.13	8.47	9.6	9.34	96.9	
M4560-04	GB8D0	4	1.17	8.47	9.64	9.59	99.4	
M4560-05	GB8D0D	5	1.17	8.47	9.64	9.59	99.4	
M4560-06	GB8D0S	6	1.17	8.47	9.64	9.59	99.4	
M4560-07	GB8D1	7	1.12	8.68	9.8	9.41	95.5	
M4560-08	GB8D2	8	1.17	8.55	9.72	9.51	97.5	
M4560-09	GB8D3	9	1.19	8.62	9.81	9.48	96.2	
M4560-10	GB8D4	10	1.18	8.45	9.63	9.25	95.5	
M4560-11	GB8D5	11	1.18	8.57	9.75	9.5	97.1	
M4560-12	GB8D6	12	1.16	8.61	9.77	9.55	97.4	
M4560-13	GB8D7	13	1.17	8.42	9.59	9.49	98.8	
M4560-14	GB8D8	14	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$