



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 11/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:20
In Date: 11/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 11/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117293

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4512-01	A4263	1	1.15	8.81	9.96	9.13	90.6	
M4512-02	A4264	2	1.13	8.71	9.84	8.85	88.6	
M4512-03	A4265	3	1.16	8.58	9.74	7.38	72.5	
M4512-04	A4266	4	1.15	8.50	9.65	8.6	87.6	
M4512-05	A4266MS	5	1.15	8.50	9.65	8.6	87.6	
M4512-06	A4266MSD	6	1.15	8.50	9.65	8.6	87.6	
M4512-07	A4267	7	1.13	8.65	9.78	8.41	84.2	
M4512-08	A4268	8	1.16	8.44	9.6	7.65	76.9	
M4512-09	A4277	9	1.14	8.49	9.63	7.4	73.7	
M4512-10	A4285	10	1.17	8.80	9.97	8.4	82.2	
M4512-11	A4286	11	1.13	8.45	9.58	8.66	89.1	
M4512-12	A4287	12	1.16	8.82	9.98	9.51	94.7	
M4512-13	A4296	13	1.15	8.81	9.96	7.87	76.3	
M4512-14	A42A0	14	1.12	8.56	9.68	7.36	72.9	
M4512-15	A42A4	15	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$