



PERCENT SOLID

Supervisor: apatel
Analyst: ketankumar
Date: 11/19/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 11/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 108
Time OUT: 09:50
Out Date: 11/19/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117436

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4702-01	DBLN1	1	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE
M4702-02	DBLN9	2	1.14	8.48	9.62	7.63	76.5	
M4702-03	DBLN9MS	3	1.14	8.48	9.62	7.63	76.5	
M4702-04	DBLN9MSD	4	1.14	8.48	9.62	7.63	76.5	
M4702-05	DBLP0	5	1.15	8.82	9.97	7.94	77.0	
M4702-06	DBLP1	6	1.17	8.45	9.62	8.26	83.9	
M4702-07	DBLP2	7	1.18	8.63	9.81	8.43	84.0	
M4702-08	DBLP3	8	1.14	8.37	9.51	7.82	79.8	
M4702-09	DBLP4	9	1.18	8.48	9.66	7.53	74.9	
M4702-10	DBLP5	10	1.16	8.52	9.68	8.62	87.6	
M4702-11	DBLP6	11	1.17	8.49	9.66	8.49	86.2	
M4702-12	DBLP7	12	1.15	8.81	9.96	8.16	79.6	
M4702-13	DBLP8	13	1.17	8.53	9.7	8.01	80.2	
M4702-14	DBLP9	14	1.17	8.44	9.61	8.32	84.7	
M4702-15	DBLQ0	15	1.19	8.67	9.86	7.44	72.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$