



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 12/2/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 12/01/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:29
Out Date: 12/02/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117601

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4872-01	DT-CON-PAD	1	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4873-01	A2-2-6.5	2	1.15	8.71	9.86	8.92	89.2	
M4873-02	A2-3-6.5	3	1.15	8.51	9.66	8.12	81.9	
M4873-03	A2-4-6.5	4	1.16	8.51	9.67	8.37	84.7	
M4873-04	A2-5-6.5	5	1.13	8.74	9.87	8.74	87.1	
M4873-05	A2-6-6.5	6	1.18	8.52	9.7	8.49	85.8	
M4873-06	A2-7-6.5	7	1.13	8.39	9.52	7.89	80.6	
M4873-07	A2-8-6.5	8	1.15	8.79	9.94	8.81	87.1	
M4873-08	A2-9-6.5	9	1.12	8.70	9.82	8.56	85.5	
M4873-09	A6-16-7.5	10	1.19	8.69	9.88	8.62	85.5	
M4873-10	A6-17-7.5	11	1.18	8.48	9.66	7.68	76.7	
M4873-11	A6-18-4.5	12	1.13	8.45	9.58	8.7	89.6	
M4873-12	A6-19-4.5	13	1.14	8.81	9.95	9.3	92.6	
M4873-13	A6-20-7.5	14	1.18	8.79	9.97	8.27	80.7	
M4876-01	CL-02-120121	15	1.14	8.51	9.65	9.14	94.0	
M4876-02	CL-02-120121-EPH-2	16	1.11	8.58	9.69	8.92	91.0	
M4877-01	PL-2-120121	17	1.16	8.80	9.96	8.89	87.8	
M4877-02	PL-2-120121-EPH-2	18	1.14	8.79	9.93	8.95	88.9	
M4877-03	PL-5-120121	19	1.17	8.60	9.77	8.88	89.7	
M4877-04	PL-5-120121-EPH-2	20	1.18	8.75	9.93	9.18	91.4	
M4891-01	MINEBROOK-1-COMP	21	1.11	8.77	9.88	8.77	87.3	
M4891-02	MINEBROOK-1-EPH-2	22	1.13	8.77	9.9	8.15	80.0	
M4891-03	MINEBROOK-2-COMP	23	1.11	8.50	9.61	8.46	86.5	
M4891-04	MINEBROOK-2-EPH-2	24	1.16	8.38	9.54	8.28	85.0	
M4893-01	EO-1-120121	25	1.11	8.41	9.52	8.67	89.9	
M4893-02	EO-1-120121-EPH-2	26	1.18	8.39	9.57	8.77	90.5	
M4893-02D	EO-1-120121-EPH-2DUP	27	1.19	8.41	9.6	8.79	90.4	
M4894-01	T-5G	28	1.14	8.82	9.96	9.2	91.4	



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 12/2/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 12/01/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:29
Out Date: 12/02/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117601

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4895-01	SLAB-1	29	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-02	SLAB-2	30	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-03	SLAB-3	31	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-04	SLAB-4	32	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-05	SLAB-5	33	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-06	SLAB-6	34	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-07	SLAB-7	35	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-08	SLAB-8	36	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-09	SLAB-9	37	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-10	SLAB-10	38	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M4895-11	SLAB-11	39	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$