



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 12/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 12/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117694

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4962-01	KMH293R-1-1	1	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-02	KMH293R-1-2	2	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-03	HIA775H-1-1	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-04	HIA775H-1-2	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-05	HIA775H-2-1	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-06	HIA775H-2-2	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-07	MOO-2021-26-A	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-08	MOO-2021-26-B	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-09	MOO-2021-27-A	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-10	MOO-2021-27-B	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-11	MOO-2021-28-A	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-12	MOO-2021-28-B	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-13	MOO-2021-29-A	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-14	MOO-2021-29-B	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-15	MOO-2021-30-A	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-16	MOO-2021-31-A	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-17	MOO-2021-31-B	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-18	MOO-2021-32-A	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-19	MOO-2021-32-B	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-20	MOO-2021-33-A	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-21	MOO-2021-34-A	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-22	MOO-2021-35-A	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-23	MOO-2021-36-A	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-24	MOO-2021-37-A	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-25	MOO-2021-37-B	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-26	MOO-2021-38-A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-27	MOO-2021-38-B	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4962-28	MOO-2021-39-A	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 12/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 12/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117694

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4962-29	MOO-2021-39-B	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
M4966-01	P5-2	30	1.15	8.77	9.92	8.89	88.3	
M4966-03	P5-10	31	1.13	8.77	9.9	8.3	81.8	
M4966-04	P4-11	32	1.16	8.64	9.8	8.55	85.5	
M4966-05	P1-COMP	33	1.17	8.73	9.9	8.48	83.7	
M4966-07	P1-8	34	1.14	8.37	9.51	8.68	90.1	
M4966-08	P6-7	35	1.18	8.74	9.92	9.08	90.4	
M4966-09	P7-4	36	1.16	8.54	9.7	8.27	83.3	
M4966-10	P8-8	37	1.12	8.68	9.8	8.72	87.6	
M4966-11	P9-3	38	1.17	8.73	9.9	8.46	83.5	
M4966-12	P10-11	39	1.15	8.81	9.96	8.37	82.0	
M4966-14	P2-COMP	40	1.14	8.79	9.93	8.76	86.7	
M4966-16	P12-4	41	1.13	8.46	9.59	8.45	86.5	
M4966-17	P13-6	42	1.15	8.82	9.97	8.96	88.5	
M4966-18	P14-8	43	1.14	8.78	9.92	8.51	83.9	
M4966-19	P14-11	44	1.19	8.68	9.87	8.77	87.3	
M4966-21	P15-7	45	1.16	8.59	9.75	8.32	83.4	
M4966-22	P3-COMP	46	1.16	8.81	9.97	8.77	86.4	
M4966-24	P16-5	47	1.18	8.72	9.9	8.4	82.8	
M4966-25	P18-4	48	1.18	8.63	9.81	8.56	85.5	
M4966-26	P19-8	49	1.15	8.82	9.97	9.44	94.0	
M4966-27	P20-8	50	1.18	8.78	9.96	8.54	83.8	
M4966-29	P4-COMP	51	1.16	8.61	9.77	8.52	85.5	
M4967-01	TP-7S	52	1.13	8.47	9.6	8.19	83.4	
M4967-03	TP-7S-EPH	53	1.14	8.80	9.94	8.68	85.7	
M4967-04	TP-7D	54	1.15	8.82	9.97	8.98	88.8	
M4967-06	TP-7D-EPH	55	1.15	8.82	9.97	9.38	93.3	
M4967-07	TP-8S	56	1.11	8.47	9.58	8.41	86.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 12/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 12/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117694

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M4967-09	TP-8S-EPH	57	1.11	8.66	9.77	8.65	87.1	
M4967-10	TP-11S	58	1.19	8.58	9.77	8.63	86.7	
M4967-12	TP-11S-EPH	59	1.18	8.79	9.97	8.95	88.4	
M4967-13	TP-11D	60	1.11	8.58	9.69	8.79	89.5	
M4967-15	TP-11D-EPH	61	1.13	8.63	9.76	8.84	89.3	
M4967-15D	TP-11D-EPHDUP	62	1.14	8.65	9.79	8.86	89.2	
M4973-01	HD-02-120821	63	1.15	8.80	9.95	9.49	94.8	
M4973-02	HD-02-120821-EPH-2	64	1.19	8.49	9.68	9.24	94.8	
M4980-01	OK-03-120821	65	1.15	8.80	9.95	9.63	96.4	
M4980-02	OK-03-120821-EPH-2	66	1.17	8.63	9.8	9.33	94.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$