



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:15
In Date: 12/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117724

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5000-01	A0BA4	1	1.14	8.63	9.77	9.04	91.5	
M5000-02	A0BA4D	2	1.14	8.63	9.77	9.04	91.5	
M5000-03	A0BA4S	3	1.14	8.63	9.77	9.04	91.5	
M5000-04	A0BA5	4	1.14	8.49	9.63	8.94	91.9	
M5000-05	A0BA6	5	1.17	8.78	9.95	9.34	93.1	
M5000-06	A0BA7	6	1.12	8.45	9.57	8.67	89.3	
M5000-07	A0BA8	7	1.16	8.81	9.97	9.6	95.8	
M5000-08	A0BA9	8	1.15	8.83	9.98	9.29	92.2	
M5000-09	A0BB0	9	1.12	8.62	9.74	8.43	84.8	
M5000-10	A0BB1	10	1.18	8.45	9.63	8.83	90.5	
M5000-11	A0BB2	11	1.19	8.58	9.77	9.06	91.7	
M5000-12	A0BB3	12	1.18	8.59	9.77	9.2	93.4	
M5000-13	A0BC4	13	1.12	8.63	9.75	8.63	87.0	
M5000-14	A0BC5	14	1.15	8.39	9.54	8.53	88.0	
M5000-15	A0BC6	15	1.17	8.75	9.92	9.34	93.4	
M5000-16	A0BC7	16	1.19	8.46	9.65	7.02	68.9	
M5000-17	A0BC8	17	1.13	8.44	9.57	8.55	87.9	
M5000-19	A0BD1	18	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$