



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:35
In Date: 12/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117726

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5004-01	BGKW0	1	1.11	8.66	9.77	9.05	91.7	
M5004-02	BGKW1	2	1.16	8.51	9.67	8.58	87.2	
M5004-03	BGKW2	3	1.15	8.83	9.98	8.72	85.7	
M5004-04	BGKW3	4	1.17	8.81	9.98	8.3	80.9	
M5004-05	BGKW4	5	1.11	8.67	9.78	7.34	71.9	
M5004-06	BGKW5	6	1.13	8.45	9.58	6.95	68.9	
M5004-07	BGKW6	7	1.16	8.49	9.65	7.23	71.5	
M5004-08	BGKW7	8	1.15	8.82	9.97	9.4	93.5	
M5004-09	BGKW8	9	1.15	8.82	9.97	9.08	89.9	
M5004-10	BGKW9	10	1.16	8.82	9.98	8.85	87.2	
M5004-11	BGKX0	11	1.13	8.64	9.77	7.32	71.6	
M5004-12	BGKX1	12	1.17	8.68	9.85	7.4	71.8	
M5004-14	BGKX3	13	1.13	8.54	9.67	8.33	84.3	
M5004-15	BGKX4	14	1.15	8.67	9.82	8.66	86.6	
M5004-16	BGKX5	15	1.16	8.82	9.98	8.73	85.8	
M5004-17	BGKX6	16	1.14	8.84	9.98	7.47	71.6	
M5004-18	VHBLK001	17	1.00	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$