



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:26
Out Date: 12/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5017-01	TP1-E1	1	1.14	8.79	9.93	9.09	90.4	
M5017-02	TP1-E2	2	1.18	8.43	9.61	9.13	94.3	
M5017-03	TP1-E3	3	1.17	8.38	9.55	9.05	94.0	
M5017-04	TP1-E4	4	1.17	8.80	9.97	9.42	93.7	
M5017-05	TP1-E5	5	1.15	8.82	9.97	8.37	81.9	
M5017-06	TP1-E6	6	1.14	8.80	9.94	8.43	82.8	
M5017-07	TP1-E7	7	1.18	8.79	9.97	8.48	83.0	
M5017-08	TP1-E8	8	1.14	8.74	9.88	8.12	79.9	
M5017-09	TP1	9	1.16	8.81	9.97	9.04	89.4	
M5017-11	TP2-E1	10	1.17	8.79	9.96	9.42	93.9	
M5017-12	TP2-E2	11	1.15	8.35	9.5	8.61	89.3	
M5017-13	TP2-E3	12	1.13	8.80	9.93	9.38	93.7	
M5017-14	TP2-E4	13	1.18	8.43	9.61	9.14	94.4	
M5017-15	TP2-E5	14	1.14	8.49	9.63	9.03	92.9	
M5017-16	TP2-E6	15	1.14	8.49	9.63	8.54	87.2	
M5017-17	TP2-E7	16	1.13	8.72	9.85	8.44	83.8	
M5017-18	TP2-E8	17	1.14	8.82	9.96	8.54	83.9	
M5017-19	TP2	18	1.19	8.65	9.84	9.02	90.5	
M5017-21	TP3-E1	19	1.13	8.61	9.74	8.95	90.8	
M5017-22	TP3-E2	20	1.19	8.78	9.97	8.82	86.9	
M5017-23	TP3-E3	21	1.16	8.65	9.81	9.05	91.2	
M5017-24	TP3-E4	22	1.19	8.71	9.9	9.05	90.2	
M5017-25	TP3-E5	23	1.13	8.38	9.51	8.62	89.4	
M5017-26	TP3-E6	24	1.19	8.53	9.72	8.54	86.2	
M5017-27	TP3-E7	25	1.15	8.40	9.55	8.22	84.2	
M5017-28	TP3-E8	26	1.19	8.53	9.72	9.24	94.4	
M5017-29	TP3	27	1.14	8.84	9.98	9.02	89.1	
M5017-31	SP1-A-E1	28	1.11	8.87	9.98	8.96	88.5	



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:26
Out Date: 12/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5017-32	SP1-A-E2	29	1.14	8.83	9.97	9.7	96.9	
M5017-33	SP1-A-E3	30	1.19	8.65	9.84	9.18	92.4	
M5017-34	SP1-A-E4	31	1.16	8.80	9.96	9.37	93.3	
M5017-35	SP1-A-E5	32	1.13	8.64	9.77	9.06	91.8	
M5017-36	SP1-A-E6	33	1.15	8.82	9.97	9.69	96.8	
M5017-37	SP1-A-E7	34	1.17	8.81	9.98	9.24	91.6	
M5017-38	SP1-A-E8	35	1.18	8.52	9.7	9.08	92.7	
M5017-39	SP1A	36	1.13	8.50	9.63	9.12	94.0	
M5017-41	SP1-B-E1	37	1.19	8.73	9.92	9.25	92.3	
M5017-42	SP1-B-E2	38	1.19	8.44	9.63	9.00	92.5	
M5017-43	SP1-B-E3	39	1.17	8.47	9.64	9.05	93.0	
M5017-44	SP1-B-E4	40	1.12	8.43	9.55	8.94	92.8	
M5017-45	SP1-B-E5	41	1.13	8.49	9.62	9.06	93.4	
M5017-46	SP1-B-E6	42	1.14	8.70	9.84	9.16	92.2	
M5017-47	SP1-B-E7	43	1.11	8.39	9.5	8.85	92.3	
M5017-48	SP1-B-E8	44	1.19	8.71	9.9	9.28	92.9	
M5017-49	SP1-B	45	1.15	8.48	9.63	9.06	93.3	
M5017-51	SP2-A-E1	46	1.19	8.68	9.87	8.9	88.8	
M5017-52	SP2-A-E2	47	1.13	8.56	9.69	8.91	90.9	
M5017-53	SP2-A-E3	48	1.18	8.77	9.95	9.22	91.7	
M5017-54	SP2-A-E4	49	1.19	8.65	9.84	9.11	91.6	
M5017-55	SP2-A-E5	50	1.16	8.82	9.98	9.15	90.6	
M5017-56	SP2-A-E6	51	1.13	8.47	9.6	8.74	89.8	
M5017-57	SP2-A-E7	52	1.13	8.37	9.5	8.68	90.2	
M5017-58	SP2-A-E8	53	1.13	8.47	9.6	8.68	89.1	
M5017-59	SP2-A	54	1.13	8.64	9.77	8.94	90.4	
M5018-01	SP2B-E1	55	1.13	8.77	9.9	8.99	89.6	
M5018-02	SP2B-E2	56	1.13	8.68	9.81	8.92	89.7	



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:26
Out Date: 12/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5018-03	SP2B-E3	57	1.13	8.81	9.94	8.83	87.4	
M5018-04	SP2B-E4	58	1.13	8.49	9.62	8.76	89.9	
M5018-05	SP2B-E5	59	1.13	8.77	9.9	9.28	92.9	
M5018-06	SP2B-E6	60	1.13	8.80	9.93	9.55	95.7	
M5018-07	SP2B-E7	61	1.13	8.45	9.58	9.03	93.5	
M5018-08	SP2B-E8	62	1.13	8.58	9.71	8.97	91.4	
M5018-09	SP2B	63	1.13	8.66	9.79	9.12	92.3	
M5018-11	TP4-E1	64	1.13	8.64	9.77	8.89	89.8	
M5018-12	TP4-E2	65	1.13	8.70	9.83	9.06	91.1	
M5018-13	TP4-E3	66	1.13	8.57	9.7	9.25	94.7	
M5018-14	TP4-E4	67	1.13	8.73	9.86	9.43	95.1	
M5018-15	TP4-E5	68	1.13	8.77	9.9	9.47	95.1	
M5018-16	TP4-E6	69	1.13	8.61	9.74	9.24	94.2	
M5018-17	TP4-E7	70	1.13	8.47	9.6	9.13	94.5	
M5018-18	TP4-E8	71	1.13	8.73	9.86	9.33	93.9	
M5018-19	TP4	72	1.13	8.69	9.82	9.27	93.7	
M5018-21	TP5-E1	73	1.13	8.69	9.82	8.95	90.0	
M5018-22	TP5-E2	74	1.13	8.84	9.97	9.66	96.5	
M5018-23	TP5-E3	75	1.13	8.75	9.88	9.53	96.0	
M5018-24	TP5-E4	76	1.13	8.68	9.81	9.41	95.4	
M5018-25	TP5-E5	77	1.13	8.77	9.9	9.72	97.9	
M5018-26	TP5-E6	78	1.13	8.69	9.82	9.38	94.9	
M5018-27	TP5-E7	79	1.13	8.85	9.98	9.69	96.7	
M5018-28	TP5-E8	80	1.13	8.76	9.89	9.47	95.2	
M5018-29	TP5	81	1.13	8.54	9.67	9.4	96.8	
M5018-31	TP6-E1	82	1.13	8.70	9.83	9.18	92.5	
M5018-32	TP6-E2	83	1.13	8.47	9.6	8.84	91.0	
M5018-33	TP6-E3	84	1.13	8.84	9.97	9.38	93.3	



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:26
Out Date: 12/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5018-34	TP6-E4	85	1.13	8.37	9.5	8.89	92.7	
M5018-35	TP6-E5	86	1.13	8.84	9.97	9.63	96.2	
M5018-36	TP6-E6	87	1.13	8.69	9.82	9.42	95.4	
M5018-37	TP6-E7	88	1.13	8.81	9.94	9.52	95.2	
M5018-38	TP6-E8	89	1.13	8.59	9.72	9.6	98.6	
M5018-39	TP6	90	1.13	8.63	9.76	9.35	95.2	
M5018-41	TP7-E1	91	1.13	8.63	9.76	8.85	89.5	
M5018-42	TP7-E2	92	1.13	8.57	9.7	9.31	95.4	
M5018-43	TP7-E3	93	1.13	8.74	9.87	9.44	95.1	
M5018-44	TP7-E4	94	1.13	8.37	9.5	9.24	96.9	
M5018-45	TP7-E5	95	1.13	8.72	9.85	9.53	96.3	
M5018-46	TP7-E6	96	1.13	8.40	9.53	9.13	95.2	
M5018-47	TP7-E7	97	1.13	8.58	9.71	9.42	96.6	
M5018-48	TP7-E8	98	1.13	8.84	9.97	9.42	93.8	
M5018-49	TP7	99	1.13	8.57	9.7	9.29	95.2	
M5018-51	TP8-E1	100	1.13	8.59	9.72	9.38	96.0	
M5018-52	TP8-E2	101	1.13	8.82	9.95	9.54	95.4	
M5018-53	TP8-E3	102	1.13	8.75	9.88	9.41	94.6	
M5018-54	TP8-E4	103	1.13	8.84	9.97	9.73	97.3	
M5018-55	TP8-E5	104	1.13	8.78	9.91	9.58	96.2	
M5018-56	TP8-E6	105	1.13	8.67	9.8	9.66	98.4	
M5018-57	TP8-E7	106	1.13	8.44	9.57	9.46	98.7	
M5018-58	TP8-E8	107	1.13	8.49	9.62	9.47	98.2	
M5018-59	TP8	108	1.13	8.59	9.72	9.44	96.7	
M5018-61	TP9-E1	109	1.13	8.52	9.65	9.22	95.0	
M5018-62	TP9-E2	110	1.13	8.72	9.85	9.5	96.0	
M5018-63	TP9-E3	111	1.13	8.84	9.97	9.65	96.4	
M5018-64	TP9-E4	112	1.13	8.37	9.5	9.14	95.7	



PERCENT SOLID

Supervisor: ketankumar

Analyst: jignesh

Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:50

In Date: 12/10/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:26

Out Date: 12/11/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5018-65	TP9-E5	113	1.13	8.42	9.55	9.03	93.8	
M5018-66	TP9-E6	114	1.13	8.70	9.83	9.54	96.7	
M5018-67	TP9-E7	115	1.13	8.82	9.95	9.36	93.3	
M5018-68	TP9-E8	116	1.13	8.46	9.59	9.24	95.9	
M5018-69	TP9	117	1.13	8.76	9.89	9.74	98.3	
M5022-01	TP10-E1	118	1.12	8.70	9.82	9.44	95.6	
M5022-02	TP10-E2	119	1.16	8.60	9.76	9.27	94.3	
M5022-03	TP10-E3	120	1.18	8.34	9.52	9.00	93.8	
M5022-04	TP10-E4	121	1.15	8.81	9.96	9.64	96.4	
M5022-05	TP10-E5	122	1.16	8.34	9.5	9.04	94.5	
M5022-06	TP10-E6	123	1.18	8.54	9.72	8.86	89.9	
M5022-07	TP10-E7	124	1.18	8.62	9.8	9.18	92.8	
M5022-08	TP10-E8	125	1.18	8.62	9.8	9.17	92.7	
M5022-09	TP10	126	1.19	8.72	9.91	9.33	93.3	
M5022-11	TP11-E1	127	1.19	8.43	9.62	9.23	95.4	
M5022-12	TP11-E2	128	1.14	8.53	9.67	9.18	94.3	
M5022-13	TP11-E3	129	1.18	8.67	9.85	9.35	94.2	
M5022-14	TP11-E4	130	1.16	8.72	9.88	9.44	95.0	
M5022-15	TP11-E5	131	1.13	8.85	9.98	9.52	94.8	
M5022-16	TP11-E6	132	1.15	8.60	9.75	9.07	92.1	
M5022-17	TP11-E7	133	1.14	8.82	9.96	9.46	94.3	
M5022-18	TP11-E8	134	1.17	8.65	9.82	9.06	91.2	
M5022-19	TP11	135	1.17	8.34	9.51	8.86	92.2	
M5022-21	TP12-E1	136	1.17	8.45	9.62	8.88	91.2	
M5022-22	TP12-E2	137	1.14	8.84	9.98	9.14	90.5	
M5022-23	TP12-E3	138	1.1	8.45	9.55	8.84	91.6	
M5022-24	TP12-E4	139	1.1	8.68	9.78	9.17	93.0	
M5022-25	TP12-E5	140	1.14	8.36	9.5	8.75	91.0	



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:50
In Date: 12/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:26
Out Date: 12/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5022-26	TP12-E6	141	1.16	8.75	9.91	8.71	86.3	
M5022-27	TP12-E7	142	1.19	8.52	9.71	8.92	90.7	
M5022-28	TP12-E8	143	1.14	8.38	9.52	8.4	86.6	
M5022-29	TP12	144	1.19	8.36	9.55	8.73	90.2	
M5022-31	TP13-E1	145	1.11	8.60	9.71	8.96	91.3	
M5022-32	TP13-E2	146	1.19	8.57	9.76	9.06	91.8	
M5022-33	TP13-E3	147	1.19	8.50	9.69	8.96	91.4	
M5022-34	TP13-E4	148	1.16	8.81	9.97	9.12	90.4	
M5022-35	TP13-E5	149	1.18	8.52	9.7	8.92	90.8	
M5022-36	TP13-E6	150	1.11	8.49	9.6	8.8	90.6	
M5022-37	TP13-E7	151	1.18	8.80	9.98	8.93	88.1	
M5022-38	TP13-E8	152	1.14	8.83	9.97	9.12	90.4	
M5022-39	TP13	153	1.15	8.51	9.66	8.89	91.0	
M5022-41	TP14-E1	154	1.15	8.40	9.55	8.87	91.9	
M5022-42	TP14-E2	155	1.15	8.80	9.95	9.22	91.7	
M5022-43	TP14-E3	156	1.13	8.67	9.8	9.09	91.8	
M5022-44	TP14-E4	157	1.17	8.46	9.63	8.82	90.4	
M5022-45	TP14-E5	158	1.19	8.51	9.7	8.89	90.5	
M5022-46	TP14-E6	159	1.13	8.45	9.58	8.86	91.5	
M5022-47	TP14-E7	160	1.1	8.79	9.89	9.11	91.1	
M5022-48	TP14-E8	161	1.19	8.63	9.82	9.19	92.7	
M5022-49	TP14	162	1.18	8.52	9.7	8.99	91.7	
M5022-51	TP15-E1	163	1.13	8.45	9.58	9.06	93.8	
M5022-52	TP15-E2	164	1.18	8.79	9.97	9.44	94.0	
M5022-53	TP15-E3	165	1.18	8.43	9.61	9.11	94.1	
M5022-54	TP15-E4	166	1.18	8.53	9.71	9.31	95.3	
M5022-55	TP15-E5	167	1.16	8.57	9.73	9.06	92.2	
M5022-56	TP15-E6	168	1.11	8.71	9.82	8.86	89.0	



PERCENT SOLID

Supervisor: ketankumar

Analyst: jignesh

Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:50

In Date: 12/10/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:26

Out Date: 12/11/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5022-57	TP15-E7	169	1.14	8.65	9.79	8.74	87.9	
M5022-58	TP15-E8	170	1.15	8.69	9.84	9.14	91.9	
M5022-59	TP15	171	1.15	8.81	9.96	9.26	92.1	
M5029-01	PSC-124035	182	1.14	8.83	9.97	6.83	64.4	GEL MATRIX
M5029-03	PSC-124019	183	1.11	8.78	9.89	8.17	80.4	
M5030-01	SU-02-121021	184	1.14	8.42	9.56	8.67	89.4	
M5030-02	SU-02-121021-EPH-2	185	1.17	8.63	9.8	8.9	89.6	
M5030-03	SU-04-121021	186	1.12	8.40	9.52	8.91	92.7	
M5030-04	SU-04-121021-EPH-2	187	1.17	8.46	9.63	8.99	92.4	
M5031-05	SVOC-GPC-BLANK	172	1.00	1.00	2.00	2.00	100.0	
M5031-06	PEST-GPC-BLANK	173	1.00	1.00	2.00	2.00	100.0	
M5031-07	PEST-GPC-BLANK-SPIKE	174	1.00	1.00	2.00	2.00	100.0	
M5031-08	PCB-GPC-BLANK	175	1.00	1.00	2.00	2.00	100.0	
M5031-09	PCB-GPC-BLANK-SPIKE	176	1.00	1.00	2.00	2.00	100.0	
M5031-10	SVOC-GPC2-BLANK	177	1.00	1.00	2.00	2.00	100.0	
M5031-11	PEST-GPC2-BLANK	178	1.00	1.00	2.00	2.00	100.0	
M5031-12	PEST-GPC2-BLANK-SPIKE	179	1.00	1.00	2.00	2.00	100.0	
M5031-13	PCB-GCP2-BLANK	180	1.00	1.00	2.00	2.00	100.0	
M5031-14	PCB-GCP2-BLANK-SPIKE	181	1.00	1.00	2.00	2.00	100.0	
M5032-01	TR-04-121021	188	1.18	8.43	9.61	9.21	95.3	
M5032-02	TR-04-121021-EPH-2	189	1.15	8.78	9.93	9.4	94.0	
M5033-01	CL-01-121021	190	1.11	8.72	9.83	9.09	91.5	
M5033-02	CL-01-121021-EPH-2	191	1.14	8.84	9.98	9.12	90.3	
M5044-01	DA-1-4.5	196	1.12	8.47	9.59	8.22	83.8	
M5044-02	DA-2-4.5	197	1.11	8.62	9.73	8.47	85.4	
M5044-03	DA-3-4.5	198	1.11	8.53	9.64	8.81	90.3	
M5044-04	DA-4-4.5	199	1.14	8.72	9.86	8.81	88.0	
M5044-05	DA-5-4.5	200	1.14	8.79	9.93	8.57	84.5	



PERCENT SOLID

Supervisor: ketankumar

Analyst: jignesh

Date: 12/11/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:50

In Date: 12/10/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:26

Out Date: 12/11/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB117738

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5044-06	DA-6-4.5	201	1.18	8.51	9.69	8.87	90.4	
M5044-07	DA-7-4.5	202	1.16	8.59	9.75	8.96	90.8	
M5044-08	DA-8-4.5	203	1.13	8.39	9.52	8.71	90.3	
M5044-09	DA-9-4.5	204	1.13	8.45	9.58	8.51	87.3	
M5044-10	DA-10-4.5	205	1.15	8.51	9.66	7.76	77.7	
M5044-11	DA-11-4.5	206	1.15	8.51	9.66	7.36	73.0	
M5044-12	DA-12-4.5	207	1.17	8.81	9.98	7.94	76.8	
M5044-13	DA-13-4.5	208	1.16	8.80	9.96	8.95	88.5	
M5044-14	DA-14-4.5	209	1.14	8.60	9.74	8.95	90.8	
M5044-15	DA-15-4.5	210	1.15	8.81	9.96	9.5	94.8	
M5044-16	DA-16-4.5	211	1.14	8.83	9.97	9.06	89.7	
M5044-17	DA-17-4.5	212	1.19	8.59	9.78	8.99	90.8	
M5044-18	DA-18-5.5	213	1.18	8.45	9.63	7.26	72.0	
M5044-19	DA-19-5.0	214	1.11	8.68	9.79	8.02	79.6	
M5044-20	DA-20-5.5	215	1.17	8.81	9.98	7.67	73.8	
M5044-21	DA-21-4.5	216	1.14	8.81	9.95	9.38	93.5	
M5051-01	P-1	192	1.14	8.60	9.74	9.06	92.1	
M5051-02	P-2	193	1.19	8.49	9.68	8.79	89.5	
M5051-03	P-3	194	1.14	8.81	9.95	9.28	92.4	
M5051-04	A5	195	1.1	8.84	9.94	8.43	82.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$