



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/14/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 12/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 12/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117770

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5064-01	SASP2-121-1-(0-5)	4	1.17	8.71	9.88	8.61	85.4	
M5064-02	SASP2-121-1-(5-10)	5	1.16	8.64	9.8	9.44	95.8	
M5064-03	SASP2-121-1-(10-17)	6	1.11	8.56	9.67	7.94	79.8	
M5064-04	SASP2-121-1-(17-32)	7	1.13	8.44	9.57	7.78	78.8	
M5064-04D	SASP2-121-1-(17-32) DUP	8	1.16	8.43	9.59	7.8	78.8	
M5065-01	UV-1WC	1	1.14	8.58	9.72	7.72	76.7	
M5067-01	BASE	2	1.15	8.38	9.53	9.27	96.9	
M5067-04	BRICK	3	1.17	8.63	9.8	9.73	99.2	
M5068-01	PL-01-121321	43	1.12	8.44	9.56	8.22	84.1	
M5068-02	PL-01-121321-EPH2	44	1.16	8.43	9.59	8.34	85.2	
M5070-01	CLI-2021-8-A	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-02	CLI-2021-9-A	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-03	CLI-2021-9-B	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-04	CLI-2021-10-A	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-05	CLI-2021-10-B	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-06	CLI-2021-11-A	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-07	CLI-2021-11-B	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-08	CLI-2021-12-A	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-09	CLI-2021-12-B	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-10	CLI-2021-13-A	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-11	CLI-2021-13-B	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-12	CLI-2021-14-A	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-13	CLI-2021-14-B	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-14	CLI-2021-15-A	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-15	CLI-2021-15-B	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-16	CLI-2021-16-A	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-17	CLI-2021-16-B	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-18	CLI-2021-17-A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 12/14/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 12/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 12/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB117770

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M5070-19	CLI-2021-17-B	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-20	CLI-2021-18-A	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-21	CLI-2021-19-A	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-22	CLI-2021-19-B	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-23	CLI-2021-20-A	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-24	CLI-2021-20-B	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-25	CLI-2021-21-A	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-26	CLI-2021-20-B	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-27	CLI-2021-22-A	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5070-28	CLI-2021-22-B	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-01	KMA666P-1-1	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-02	KMA666P-1-2	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-03	KMH430D-1-1	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-04	KMH430D-1-2	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-05	HED691G-1-1	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5071-06	HED691G-1-2	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
M5073-01	EO-02-121321	45	1.18	8.63	9.81	8.92	89.7	
M5073-02	EO-02-121321-EPH-2	46	1.15	8.61	9.76	8.91	90.1	
M5073-04	EO-03-121321	47	1.17	8.44	9.61	8.76	89.9	
M5073-05	EO-03-121321-EPH-2	48	1.12	8.74	9.86	8.99	90.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$