



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 1/3/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 01/01/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 01/02/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118038

Table with 9 columns: Lab ID, Client SampleID, Dish #, Dish Wt (g) (A), Sample Wt (g), Dish + Sample Wt (g) (B), Dish+Dry Sample Wt (g) (C), % Solid, Comments. It contains 25 rows of data for various samples.

% Solid = (C-A) \* 100 / (B-A)