



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/12/2022

OVENTEMP IN Celsius(°C): 106
Time IN: 17:15
In Date: 01/11/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:07
Out Date: 01/12/2022
Weight Check 1.0g: 0.01
Weight Check 10g: 1.00
BalanceID: M SC-4
Thermometer ID:

QC:LB118153

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1001-06	PH-2/3-3.5	1	1.12	8.40	9.52	8.51	88.0	
N1001-19	PH-5/2.5-3	2	1.15	8.81	9.96	9.47	94.4	
N1086-01	BBC	3	1.17	8.55	9.72	9.66	99.3	
N1096-01	SO-E3-02-6	4	1.13	8.70	9.83	8.87	89.0	
N1096-02	SO-E3-01-6	5	1.16	8.66	9.82	8.76	87.8	
N1096-03	SO-E7-02-6	6	1.16	8.60	9.76	7.99	79.4	
N1096-04	SO-E7-01-6	7	1.14	8.83	9.97	8.59	84.4	
N1096-05	SO-E8-01-6	8	1.11	8.62	9.73	7.56	74.8	
N1096-06	SO-E8-02-6	9	1.15	8.78	9.93	8.27	81.1	
N1096-07	SO-E2-01-8	10	1.11	8.71	9.82	8.27	82.2	
N1096-08	SO-E1-02-7	11	1.13	8.54	9.67	8.86	90.5	
N1096-08D	SO-E1-02-7DUP	12	1.15	8.54	9.69	8.87	90.4	
N1097-01	TP-1S	52	1.14	8.69	9.83	8.11	80.2	
N1097-03	TP-1S-EPH-2	53	1.17	8.63	9.8	7.65	75.1	
N1097-04	TP-1D	54	1.13	8.41	9.54	7.35	74.0	
N1097-06	TP-1D-EPH-2	55	1.15	8.81	9.96	9.1	90.2	
N1097-07	TP-2S	56	1.13	8.47	9.6	7.64	76.9	
N1097-09	TP-2S-EPH-2	57	1.17	8.51	9.68	8.36	84.5	
N1097-10	TP-2D	58	1.14	8.83	9.97	8.63	84.8	
N1097-12	TP-2D-EPH-2	59	1.15	8.79	9.94	8.42	82.7	
N1097-13	TP-3S	60	1.16	8.76	9.92	8.54	84.2	
N1097-15	TP-3S-EPH-2	61	1.18	8.76	9.94	8.96	88.8	
N1097-16	TP-3D	62	1.19	8.50	9.69	7.63	75.8	
N1097-18	TP-3D-EPH-2	63	1.14	8.83	9.97	8.86	87.4	
N1098-01	14-CON	25	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
N1099-01	VACLOT-BACKFILL-TILCON-01	13	1.14	8.76	9.9	9.35	93.7	
N1099-02	VACLOT-BACKFILL-TILCON-02	14	1.16	8.61	9.77	9.09	92.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/12/2022

OVENTEMP IN Celsius(°C): 106
Time IN: 17:15
In Date: 01/11/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:07
Out Date: 01/12/2022
Weight Check 1.0g: 0.01
Weight Check 10g: 1.00
BalanceID: M SC-4
Thermometer ID:

QC:LB118153

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1099-03	VACLOT-BACKFILL-TILCON-03	15	1.15	8.79	9.94	9.39	93.7	
N1099-04	VACLOT-BACKFILL-TILCON-04	16	1.15	8.80	9.95	9.48	94.7	
N1099-05	VACLOT-BACKFILL-TILCON-05	17	1.16	8.55	9.71	9.05	92.3	
N1099-06	VACLOT-BACKFILL-TILCON-06	18	1.16	8.76	9.92	9.27	92.6	
N1099-07	VACLOT-BACKFILL-TILCON-07	19	1.17	8.71	9.88	9.19	92.1	
N1099-08	VACLOT-BACKFILL-TILCON-08	20	1.18	8.58	9.76	9.04	91.6	
N1099-09	VACLOT-BACKFILL-TILCON-09	21	1.19	8.53	9.72	9.02	91.8	
N1099-10	VACLOT-BACKFILL-TILCON-10	22	1.15	8.40	9.55	8.86	91.8	
N1103-01	CL-01-011122	23	1.18	8.41	9.59	9.24	95.8	
N1103-02	CL-01-011122-EPH-2	24	1.14	8.51	9.65	9.36	96.6	
N1105-01	IRV-2022-19A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-02	IRV-2022-19B	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-03	IRV-2022-21A	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-04	IRV-2022-21B	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-05	IRV-2022-22A	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-06	IRV-2022-22B	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-07	IRV-2022-26A	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-08	IRV-2022-26B	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-09	IRV-2022-30A	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-10	IRV-2022-30B	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-11	IRV-2022-31A	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-12	IRV-2022-31B	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-13	IRV-2022-32A	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-14	IRV-2022-32B	39	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/12/2022

OVENTEMP IN Celsius(°C): 106
Time IN: 17:15
In Date: 01/11/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:07
Out Date: 01/12/2022
Weight Check 1.0g: 0.01
Weight Check 10g: 1.00
BalanceID: M SC-4
Thermometer ID:

QC:LB118153

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1105-15	IRV-2022-33A	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-16	IRV-2022-33B	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-17	IRV-2022-34A	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-18	IRV-2022-34B	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-19	IRV-2022-35A	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-20	IRV-2022-35B	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-21	IRV-2022-36A	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1105-22	IRV-2022-36B	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1106-01	SO-E6-02-6	64	1.14	8.77	9.91	8.67	85.9	
N1106-02	DUP-01-011122	65	1.15	8.40	9.55	8.46	87.0	
N1106-03	SO-E6-01-6	66	1.19	8.70	9.89	8.74	86.8	
N1106-04	N1106-03MS	67	1.19	8.70	9.89	8.74	86.8	
N1106-05	N1106-03MSD	68	1.19	8.70	9.89	8.74	86.8	
N1106-06	SO-0022C-01-6	69	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-07	DUP-02-011122	70	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-08	SO-0022C-02-6	71	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-09	N1106-08MS	72	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-10	N1106-08MSD	73	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-11	SO-0022C-03-6	74	1.00	1.00	2.00	2.00	100.0	encore sample receive
N1106-12	SO-A4-A41-01-6	75	1.17	8.61	9.78	8.64	86.8	
N1106-13	SO-A4-A41-02-6	76	1.11	8.47	9.58	8.8	90.8	
N1106-14	SO-A4-A41-03-6	77	1.13	8.44	9.57	8.57	88.2	
N1107-01	PT-03-011122	48	1.18	8.54	9.72	8.8	89.2	
N1107-02	PT-03-011122-EPH-2	49	1.15	8.82	9.97	9.16	90.8	
N1107-03	PT-04-011122	50	1.11	8.74	9.85	9.18	92.3	
N1107-04	PT-04-011122-EPH-2	51	1.13	8.77	9.9	8.87	88.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$