



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 01/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:45
Out Date: 01/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118319

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1217-06	72-72017	1	1.13	8.64	9.77	8.06	80.2	
N1217-07	72-72017	2	1.13	8.64	9.77	8.06	80.2	
N1217-08	VNJ-242	3	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N1217-09	VNJ-242	4	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N1229-05	SVOC-GPC-BLANK	5	1.00	1.00	2.00	2.00	100.0	
N1229-06	PEST-GPC-BLANK	6	1.00	1.00	2.00	2.00	100.0	
N1229-07	PEST-GPC-BLANK-SPIKE	7	1.00	1.00	2.00	2.00	100.0	
N1229-08	PCB-GPC-BLANK	8	1.00	1.00	2.00	2.00	100.0	
N1229-09	PCB-GPC-BLANK-SPIKE	9	1.00	1.00	2.00	2.00	100.0	
N1229-10	SVOC-GPC2-BLANK	10	1.00	1.00	2.00	2.00	100.0	
N1229-11	PEST-GPC2-BLANK	11	1.00	1.00	2.00	2.00	100.0	
N1229-12	PEST-GPC2-BLANK-SPIKE	12	1.00	1.00	2.00	2.00	100.0	
N1229-13	PCB-GCP2-BLANK	13	1.00	1.00	2.00	2.00	100.0	
N1229-14	PCB-GCP2-BLANK-SPIKE	14	1.00	1.00	2.00	2.00	100.0	
N1232-01	ELI-2022-12A	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1232-02	ELI-2022-12B	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1232-03	ELI-2022-13A	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1232-04	ELI-2022-13B	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1233-01	RT3876	42	1.14	8.84	9.98	9.09	89.9	
N1233-02	RT3876	43	1.14	8.84	9.98	9.09	89.9	
N1233-03	RT4743	44	1.19	8.62	9.81	8.39	83.5	
N1233-04	RT4743	45	1.19	8.62	9.81	8.39	83.5	
N1234-01	TR-04-01212022	46	1.15	8.81	9.96	9.00	89.1	
N1234-02	TR-04-01212022-EPH-2	47	1.11	8.77	9.88	9.18	92.0	
N1235-01	TP-1	48	1.14	8.75	9.89	8.62	85.5	
N1235-03	TP-1-EPH-2	49	1.14	8.75	9.89	8.62	85.5	
N1235-04	TP-2	50	1.19	8.72	9.91	8.44	83.1	
N1235-06	TP-2-EPH-2	51	1.19	8.72	9.91	8.44	83.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 01/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:45
Out Date: 01/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118319

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1236-01	S1-01-A	15	1.15	8.77	9.92	7.66	74.2	
N1236-04	S1-02-A	16	1.12	8.45	9.57	7.17	71.6	
N1236-07	S1-03-A	17	1.17	8.77	9.94	7.72	74.7	
N1236-09	S2-01-A	18	1.16	8.82	9.98	7.4	70.7	
N1236-12	DUP-1	19	1.11	8.84	9.95	8.24	80.7	
N1237-01	P-PARK-SCREENINGS	52	1.16	8.72	9.88	9.12	91.3	
N1237-02	P-PARK-3	53	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
N1240-01	SEC-2022-8A	20	1.00	1.00	2.00	2.00	100.0	PILC
N1240-02	SEC-2022-8B	21	1.00	1.00	2.00	2.00	100.0	PILC
N1240-03	SEC-2022-9A	22	1.00	1.00	2.00	2.00	100.0	PILC
N1240-04	SEC-2022-9B	23	1.00	1.00	2.00	2.00	100.0	PILC
N1240-05	SEC-2022-10A	24	1.00	1.00	2.00	2.00	100.0	PILC
N1240-06	SEC-2022-10B	25	1.00	1.00	2.00	2.00	100.0	PILC
N1240-07	SEC-2022-11A	26	1.00	1.00	2.00	2.00	100.0	PILC
N1240-08	SEC-2022-11B	27	1.00	1.00	2.00	2.00	100.0	PILC
N1241-01	KMA662R-1-1	28	1.00	1.00	2.00	2.00	100.0	PILC
N1241-02	KMA662R-1-2	29	1.00	1.00	2.00	2.00	100.0	PILC
N1241-03	JEB909B-1-1	30	1.00	1.00	2.00	2.00	100.0	PILC
N1241-04	JEB909B-1-2	31	1.00	1.00	2.00	2.00	100.0	PILC
N1241-05	JEC649M-1-1	32	1.00	1.00	2.00	2.00	100.0	PILC
N1241-06	JEC649M-1-2	33	1.00	1.00	2.00	2.00	100.0	PILC
N1241-07	JEC649M-2-1	34	1.00	1.00	2.00	2.00	100.0	PILC
N1241-08	JEC649M-2-2	35	1.00	1.00	2.00	2.00	100.0	PILC
N1241-09	KMA998C-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC
N1241-10	KMA998C-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC
N1242-01	PSC-124042	54	1.14	8.83	9.97	4.57	38.8	GEL MATRIX
N1242-03	PSC-124055	55	1.19	8.51	9.7	5.48	50.4	GEL MATRIX
N1242-05	PSC-124035	56	1.16	8.65	9.81	6.00	56.0	GEL MATRIX



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/24/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 18:00
In Date: 01/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:45
Out Date: 01/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118319

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1243-01	AMI-1	57	1.19	8.62	9.81	8.76	87.8	
N1243-02	AMI-1	58	1.16	8.82	9.98	8.6	84.4	
N1243-04	AMI-1-EPH-2	59	1.17	8.78	9.95	8.82	87.1	
N1243-05	AMI-2	60	1.16	8.82	9.98	8.48	83.0	
N1243-06	AMI-2	61	1.11	8.52	9.63	8.83	90.6	
N1243-08	AMI-2-EPH-2	62	1.15	8.83	9.98	8.63	84.7	
N1243-09	AMI-3	63	1.15	8.80	9.95	9.34	93.1	
N1243-10	AMI-3	64	1.11	8.51	9.62	8.58	87.8	
N1243-12	AMI-3-EPH-2	65	1.12	8.58	9.7	8.69	88.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$