



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/27/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 01/26/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:22
Out Date: 01/27/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118385

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1272-01	ETGI-348	6	1.16	8.62	9.78	9.41	95.7	
N1272-02	ETGI-348	7	1.13	8.66	9.79	8.45	84.5	
N1276-01	PCB-SF-012522-01	1	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N1276-02	PCB-SF-012522-02	2	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N1276-03	PCB-SF-012522-03	3	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N1276-04	PCB-SF-012522-04	4	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N1280-01	RMR-CON	8	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N1280-02	RMR-OILY-DRUM	9	1.00	1.00	2.00	2.00	100.0	OILY-DRUM
N1281-01	BRICK-SAMPLE	5	1.14	8.46	9.6	9.22	95.5	
N1285-01	LBS-LBS-STOCK	10	1.15	8.80	9.95	8.41	82.5	
N1285-02	LBS-LBS-STOCK	11	1.14	8.73	9.87	8.59	85.3	
N1288-01	NEB-8005-(1-1)	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1288-02	NEB-8005-(1-2)	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1288-03	NEB-8007-(1-1)	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1288-04	NEB-8007-(1-2)	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-01	FDA-693P-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-02	FDA-693P-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-03	FDA-693P-2-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-04	FDA-693P-2-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-05	CLI-2022-9A	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1289-06	CLI-2022-9B	21	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$