



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/7/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 15:50
In Date: 02/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 02/05/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118536

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1409-01	A5720	1	1.09	8.89	9.98	8.93	88.2	
N1409-02	A5721	2	1.1	8.70	9.8	9.29	94.1	
N1409-03	A5722	3	1.11	8.74	9.85	8.9	89.1	
N1409-04	A5722D	4	1.11	8.74	9.85	8.9	89.1	
N1409-05	A5722S	5	1.11	8.74	9.85	8.9	89.1	
N1409-06	A5723	6	1.11	8.74	9.85	8.83	88.3	
N1409-07	A5724	7	1.1	8.76	9.86	8.33	82.5	
N1409-08	A5725	8	1.07	8.79	9.86	7.71	75.5	
N1409-09	A5726	9	1.08	8.88	9.96	8.8	86.9	
N1409-10	A5727	10	1.12	8.75	9.87	8.39	83.1	
N1409-11	A5728	11	1.08	8.89	9.97	9.55	95.3	
N1409-12	A5729	12	1.07	8.91	9.98	8.49	83.3	
N1409-13	A5732	13	1.00	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$