



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/21/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 02/18/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 02/19/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118737

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1598-01	GS-1	1	1.18	8.63	9.81	8.99	90.5	
N1598-02	GS-2	2	1.15	8.48	9.63	8.94	91.9	
N1598-03	GS-3	3	1.15	8.60	9.75	8.96	90.8	
N1598-04	GS-4	4	1.14	8.50	9.64	8.81	90.2	
N1598-05	GS-5	5	1.19	8.58	9.77	8.92	90.1	
N1598-06	GS-6	6	1.11	8.76	9.87	9.06	90.8	
N1598-07	GS-7	7	1.15	8.82	9.97	9.04	89.5	
N1598-08	GS-8	8	1.13	8.56	9.69	9.09	93.0	
N1600-01	C0AA0	9	1.18	8.48	9.66	8.38	84.9	
N1600-02	C0AA1	10	1.15	8.38	9.53	7.86	80.1	
N1600-03	C0AA2	11	1.14	8.67	9.81	8.24	81.9	
N1600-04	C0AA4	12	1.19	8.62	9.81	8.48	84.6	
N1600-05	C0AA6	13	1.19	8.63	9.82	9.09	91.5	
N1600-06	C0AA7	14	1.15	8.59	9.74	7.95	79.2	
N1600-07	C0AA8	15	1.19	8.76	9.95	8.14	79.3	
N1600-08	C0AA9	16	1.16	8.72	9.88	8.28	81.7	
N1601-01	K22102Q-1-1	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-02	K22102Q-1-2	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-03	BC157517-1-1	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-04	BC157517-1-2	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-05	BC157517-1-3	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-06	BC157517-1-4	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-07	CLI-2022-10A	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-08	CLI-2022-11A	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-09	CLI-2022-12A	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-10	CLI-2022-13A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-11	CLI-2022-14A	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-12	CLI-2022-14B	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/21/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 02/18/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 02/19/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118737

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1601-13	CLI-2022-15A	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-14	CLI-2022-15B	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-15	CLI-2022-16A	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-16	CLI-2022-16B	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-17	CLI-2022-17A	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-18	CLI-2022-17B	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-19	CLI-2022-18A	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-20	CLI-2022-18B	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-21	CLI-2022-19A	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-22	CLI-2022-19B	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-23	CLI-2022-20A	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-24	CLI-2022-21A	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1601-25	CLI-2022-22A	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1602-01	20220075-BERKELEY-2-COMP	42	1.15	8.81	9.96	8.77	86.5	sediment sample
N1602-02	20220076-BERKELEY-3-COMP	43	1.15	8.81	9.96	8.12	79.1	sediment sample
N1602-03	20220077A-BERKELEY-4A-COMP	44	1.13	8.73	9.86	4.68	40.7	sediment sample
N1602-04	20220077B-BERKELEY-4B-COMP	45	1.15	8.78	9.93	7.4	71.2	sediment sample
N1602-05	20220078-BERKELEY-6-COMP	46	1.16	8.61	9.77	5.84	54.4	sediment sample
N1602-06	20220082-AMENDED-BERKELEY-2-COMP	47	1.18	8.55	9.73	8.22	82.3	sediment sample
N1602-07	20220083-AMENDED-BERKELEY-3-COMP	48	1.19	8.73	9.92	8.18	80.1	sediment sample
N1602-08	20220084A-AMENDED-BERKELEY-8A-COMP	49	1.16	8.65	9.81	5.15	46.1	sediment sample
N1602-09	20220084B-AMENDED-BERKELEY-4B-COMP	50	1.15	8.80	9.95	7.12	67.8	sediment sample
N1602-10	20220086-AMENDED-BERKELEY-6-COMP	51	1.19	8.50	9.69	5.92	55.6	sediment sample



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/21/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 02/18/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 02/19/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118737

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1609-01	PCB-021722-SR-01	52	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N1610-01	TR-05-021822	53	1.19	8.41	9.6	8.47	86.6	
N1610-02	TR-05-021822-EPH-2	54	1.16	8.72	9.88	8.64	85.8	
N1611-01	1746-SOIL	57	1.12	8.49	9.61	3.09	23.2	sediment sample
N1611-02	1747	58	1.12	8.47	9.59	2.52	16.5	sediment sample
N1611-03	1748	59	1.19	8.61	9.8	3.38	25.4	sediment sample
N1611-04	SOIL-COMP	60	1.19	8.79	9.98	3.43	25.5	sediment sample
N1612-01	WASTE-CHARACTER	55	1.18	8.46	9.64	8.56	87.2	
N1614-01	WASTE-BUILDING	56	1.19	8.61	9.8	8.92	89.8	
N1616-05	SVOC-GPC-BLANK	61	1.00	1.00	2.00	2.00	100.0	
N1616-06	PEST-GPC-BLANK	62	1.00	1.00	2.00	2.00	100.0	
N1616-07	PEST-GPC-BLANK-SPIKE	63	1.00	1.00	2.00	2.00	100.0	
N1616-08	PCB-GPC-BLANK	64	1.00	1.00	2.00	2.00	100.0	
N1616-09	PCB-GPC-BLANK-SPIKE	65	1.00	1.00	2.00	2.00	100.0	
N1616-10	SVOC-GPC2-BLANK	66	1.00	1.00	2.00	2.00	100.0	
N1616-11	PEST-GPC2-BLANK	67	1.00	1.00	2.00	2.00	100.0	
N1616-12	PEST-GPC2-BLANK-SPIKE	68	1.00	1.00	2.00	2.00	100.0	
N1616-13	PCB-GCP2-BLANK	69	1.00	1.00	2.00	2.00	100.0	
N1616-14	PCB-GCP2-BLANK-SPIKE	70	1.00	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$