



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/28/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/25/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 02/26/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118856

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1655-01	REG-PUMP	1	1.12	4.74	5.86	5.38	89.9	
N1655-02	REG-MID	2	1.13	4.02	5.15	4.56	85.3	
N1655-03	REG-DISP	3	1.14	4.08	5.22	4.62	85.3	
N1655-04	PREM-PUMP	4	1.19	4.09	5.28	3.86	65.3	
N1655-05	PREM-MID	5	1.15	4.89	6.04	5.26	84.0	
N1655-06	PREM-REAR	6	1.12	5.02	6.14	5.36	84.5	
N1655-07	REMOTE-1	7	1.16	6.61	7.77	6.95	87.6	
N1655-08	UNK-MID	8	1.18	6.48	7.66	6.48	81.8	
N1655-09	UNK-SOUTH	9	1.19	6.14	7.33	6.07	79.5	
N1655-10	UNK-NORTH	10	1.18	6.06	7.24	6.07	80.7	
N1664-30	TP-10-EPH2	11	1.15	8.83	9.98	9.08	89.8	
N1674-01	36CI-A	22	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N1674-02	36CI-B	23	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N1676-05	SVOC-GPC-BLANK	12	1.00	1.00	2.00	2.00	100.0	
N1676-06	PEST-GPC-BLANK	13	1.00	1.00	2.00	2.00	100.0	
N1676-07	PEST-GPC-BLANK-SPIKE	14	1.00	1.00	2.00	2.00	100.0	
N1676-08	PCB-GPC-BLANK	15	1.00	1.00	2.00	2.00	100.0	
N1676-09	PCB-GPC-BLANK-SPIKE	16	1.00	1.00	2.00	2.00	100.0	
N1676-10	SVOC-GPC2-BLANK	17	1.00	1.00	2.00	2.00	100.0	
N1676-11	PEST-GPC2-BLANK	18	1.00	1.00	2.00	2.00	100.0	
N1676-12	PEST-GPC2-BLANK-SPIKE	19	1.00	1.00	2.00	2.00	100.0	
N1676-13	PCB-GCP2-BLANK	20	1.00	1.00	2.00	2.00	100.0	
N1676-14	PCB-GCP2-BLANK-SPIKE	21	1.00	1.00	2.00	2.00	100.0	
N1678-01	OK-02-022522	24	1.17	8.52	9.69	8.26	83.2	
N1678-02	OK-02-022522-EPH-2	25	1.14	8.83	9.97	8.45	82.8	
N1679-01	RO-72-11944	26	1.18	8.63	9.81	8.43	84.0	
N1679-03	RO-72-11944-EPH-2	27	1.19	8.62	9.81	8.43	84.0	
N1679-04	RO-VNJ-221	28	1.11	8.41	9.52	7.93	81.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 2/28/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 02/25/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 02/26/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118856

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1679-06	RO-VNJ-221-EPH-2	29	1.11	8.41	9.52	7.93	81.1	
N1681-01	TP-1-1-8FT	30	1.15	8.47	9.62	8.5	86.8	
N1681-02	TP-1-1-8FT-EPH-2	31	1.15	8.47	9.62	8.5	86.8	
N1681-03	TP-1-2FT	32	1.16	8.46	9.62	8.49	86.6	
N1681-04	TP-2-1-8FT	33	1.18	8.70	9.88	8.49	84.0	
N1681-05	TP-2-1-8FT-EPH-2	34	1.18	8.70	9.88	8.49	84.0	
N1681-06	TP-2-2FT	35	1.17	8.59	9.76	8.37	83.8	
N1681-07	TP-3-1-8FT	36	1.14	8.82	9.96	8.37	82.0	
N1681-08	TP-3-1-8FT-EPH-2	37	1.14	8.82	9.96	8.37	82.0	
N1681-09	TP-3-2FT	38	1.19	8.71	9.9	8.64	85.5	
N1687-01	COMP-12-16-2021-DRUM	39	1.00	1.00	2.00	2.00	100.0	STONE-DEBRIS SAMPLE
N1687-02	11-16-21-DRUM	40	1.00	1.00	2.00	2.00	100.0	DEBRIS
N1687-03	12-29-21-DRUM	41	1.18	8.80	9.98	8.25	80.3	
N1687-04	2-17-22-DRUM	42	1.00	1.00	2.00	2.00	100.0	STONE SAMPLE, 100% SOLIDS

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$