



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/4/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/03/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1730-01	1-2-SAMPLE-LOCATION-1	1	1.15	8.80	9.95	8.97	88.9	
N1730-02	3-4-SAMPLE-LOCATION-2	2	1.14	8.79	9.93	9.05	90.0	
N1730-03	5-6-SAMPLE-LOCATION-3	3	1.12	8.44	9.56	8.63	89.0	
N1730-04	7-8-SAMPLE-LOCATION-4	4	1.19	8.58	9.77	8.79	88.6	
N1730-05	9-10-SAMPLE-LOCATION-5	5	1.19	8.78	9.97	9.17	90.9	
N1730-06	11-12-SAMPLE-LOCATION-6	6	1.19	8.63	9.82	8.72	87.3	
N1730-07	13-14-SAMPLE-LOCATION-7	7	1.14	8.41	9.55	8.51	87.6	
N1730-08	15-16-SAMPLE-LOCATION-8	8	1.11	8.55	9.66	8.48	86.2	
N1730-09	17-18-SAMPLE-LOCATION-9	9	1.14	8.70	9.84	8.71	87.0	
N1730-10	19-20-SAMPLE-LOCATION-10	10	1.15	8.81	9.96	9.16	90.9	
N1730-11	21-22-SAMPLE-LOCATION-11	11	1.19	8.61	9.8	8.79	88.3	
N1730-12	23-24-SAMPLE-LOCATION-12	12	1.19	8.69	9.88	8.63	85.6	
N1730-13	25-26-SAMPLE-LOCATION-13	13	1.16	8.81	9.97	9.12	90.4	
N1730-14	27-28-SAMPLE-LOCATION-14	14	1.16	8.72	9.88	8.79	87.5	
N1730-15	29-30-SAMPLE-LOCATION-15	15	1.19	8.76	9.95	8.85	87.4	
N1730-16	31-32-SAMPLE-LOCATION-16	16	1.15	8.40	9.55	8.58	88.5	
N1759-01	DISP-1	17	1.14	7.06	8.2	6.85	80.9	
N1759-02	DISP-2	18	1.18	6.20	7.38	5.96	77.1	
N1759-03	DISP-3	19	1.18	5.50	6.68	5.36	76.0	
N1759-04	DISP-4	20	1.19	5.59	6.78	5.88	83.9	
N1759-05	LINE-T-1	21	1.13	6.43	7.56	6.37	81.5	
N1759-06	LINE-1-2	22	1.18	6.88	8.06	6.72	80.5	
N1759-07	LINE-2-3	23	1.18	6.81	7.99	6.66	80.5	
N1759-08	LINE-3-4	24	1.12	6.27	7.39	6.27	82.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/4/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/03/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1760-01	LB-26UP-(0-0.5)-030222	25	1.14	8.82	9.96	7.5	72.1	
N1760-02	LB-26PIPE-(0-0.5)-030222	26	1.12	8.55	9.67	5.97	56.7	
N1760-04	LB-26DN(0-0.5)-030222	27	1.11	8.86	9.97	7.46	71.7	
N1760-05	DUP-1-030222	28	1.15	8.80	9.95	8.14	79.4	
N1760-06	LB-26A-(2-2.5)-030222	29	1.15	8.45	9.6	7.65	76.9	
N1760-06D	LB-26A-(2-2.5)-030222DUP	30	1.18	8.45	9.63	7.67	76.8	
N1762-01	282	31	1.17	8.61	9.78	8.54	85.6	
N1762-02	282	32	1.15	8.80	9.95	8.84	87.4	
N1774-01	DRUM-COMP-A	43	1.15	8.80	9.95	9.19	91.4	
N1774-02	DRUM-COMP-A-EPH-2	44	1.16	8.49	9.65	9.07	93.2	
N1774-03	DRUM-COMP-B	45	1.18	8.58	9.76	9.2	93.5	
N1774-04	DRUM-COMP-B-EPH-2	46	1.16	8.74	9.9	9.25	92.6	
N1785-01	TR-4-030322	47	1.11	8.69	9.8	9.19	93.0	
N1785-02	TR-4-030322-EPH-2	48	1.16	8.45	9.61	8.74	89.7	
N1786-01	HEF4380-1-2	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1786-02	HEF4380-2-2	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1789-01	CHRT-27165	51	1.14	8.81	9.95	6.43	60.0	
N1789-02	CHRT-27165-EPH-2	52	1.19	8.42	9.61	6.62	64.5	
N1790-01	DUM-17	33	1.17	8.80	9.97	8.9	87.8	
N1790-02	DUM-17-EPH-2	34	1.14	8.81	9.95	8.87	87.7	
N1790-03	DUM-18	35	1.17	8.46	9.63	8.48	86.4	
N1790-04	DUM-18-EPH-2	36	1.1	8.79	9.89	8.85	88.2	
N1790-05	DUM-19	37	1.17	8.71	9.88	8.8	87.6	
N1790-06	DUM-19-EPH-2	38	1.11	8.52	9.63	8.29	84.3	
N1790-07	DUM-20	39	1.16	8.80	9.96	9.5	94.8	
N1790-08	DUM-20-EPH-2	40	1.15	8.81	9.96	8.78	86.6	
N1790-09	DUM-21	41	1.15	8.81	9.96	8.88	87.7	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/4/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 03/03/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB118970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1790-10	DUM-21-EPH-2	42	1.18	8.75	9.93	9.43	94.3	
N1792-01	BC260952L02-1-1	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1792-02	KZ71910-1-2	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
N1792-03	KZ71910-2-2	55	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$