



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 3/29/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 03/15/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 03/16/2022
Weight Check 1.0g: 0.01
Weight Check 10g: 0.1
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119161

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1943-01	B1	1	1.15	8.81	9.96	9.00	89.1	
N1943-02	B2	2	1.15	8.78	9.93	8.26	81.0	
N1943-03	TA1	3	1.19	8.40	9.59	8.47	86.7	
N1943-04	TB1	4	1.15	8.78	9.93	8.58	84.6	
N1943-05	TB2	5	1.18	8.75	9.93	8.59	84.7	
N1943-06	DF1	6	1.14	8.77	9.91	8.62	85.3	
N1945-01	7PH-3-13.5-14	7	1.14	8.74	9.88	8.88	88.6	
N1945-02	7PH-4-13.5-14	8	1.19	8.62	9.81	8.74	87.6	
N1945-03	7PH-5-14-14.5	9	1.16	8.57	9.73	8.69	87.9	
N1945-04	7SB-1-1.5-2	10	1.18	8.56	9.74	9.14	93.0	
N1945-05	7SB-1-3-3.5	11	1.19	8.35	9.54	8.87	92.0	
N1945-06	10PH-1-11.5-12	12	1.18	8.56	9.74	9.24	94.2	
N1945-07	10PH-2-11.5-12	13	1.15	8.40	9.55	8.94	92.7	
N1945-08	4PH-1-7.5-8	14	1.11	8.49	9.6	8.56	87.8	
N1945-09	4PH-2-7.5-8	15	1.16	8.40	9.56	8.51	87.5	
N1945-10	4PH-3-7.5-8	16	1.14	8.56	9.7	8.51	86.1	
N1949-01	CP-BH2-031422-2-3	17	1.18	8.60	9.78	8.64	86.7	
N1949-02	CP-BH2-031422-5-6	18	1.18	8.46	9.64	8.91	91.4	
N1949-03	CP-BH1-031422-0-1	19	1.18	8.41	9.59	8.86	91.3	
N1949-04	CP-BH1-031422-2-3	20	1.16	8.62	9.78	9.08	91.9	
N1949-05	CP-BH1-031422-5-6	21	1.18	8.75	9.93	9.47	94.7	
N1949-06	CP-BH1-031422-5-6-DUP	22	1.18	8.47	9.65	9.28	95.6	
N1949-09	CP-BH2-031422-0-2	23	1.19	8.40	9.59	8.75	90.0	
N1949-10	CP-BH2-031422-2-4	24	1.19	8.56	9.75	9.04	91.7	
N1949-11	CP-BH2-031422-4-6	25	1.11	8.63	9.74	9.2	93.7	
N1949-12	CP-BH1-031422-0-2	26	1.18	8.78	9.96	9.13	90.5	
N1949-13	CP-BH1-031422-2-4	27	1.12	8.70	9.82	9.38	94.9	
N1949-14	CP-BH1-031422-4-6	28	1.19	8.59	9.78	9.35	95.0	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 3/29/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:45
In Date: 03/15/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 03/16/2022
Weight Check 1.0g: 0.01
Weight Check 10g: 0.1
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119161

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1949-15	CP-BH1-031422-4-6-DUP	29	1.17	8.39	9.56	9.08	94.3	
N1949-16	CP-BH1-031422-0-0.5	30	1.18	8.75	9.93	7.79	75.5	
N1949-17	CP-BH1-031422-0-0.5-DUP	31	1.19	8.56	9.75	8.45	84.8	
N1949-18	CP-BH2-031422-0-1	32	1.17	8.44	9.61	8.56	87.6	
N1951-01	DB-27-(22-25)-3-12-22	33	1.18	8.63	9.81	7.68	75.3	
N1951-02	DB-27-(15-20)-3-12-22	34	1.17	8.48	9.65	7.46	74.2	
N1951-03	DB-15-(10-15)-3-12-22	35	1.16	8.80	9.96	7.87	76.3	
N1951-04	DB-15-(20-25)-3-12-22	36	1.18	8.53	9.71	7.52	74.3	
N1951-05	DB-14-(13-15)-3-12-22	37	1.19	8.70	9.89	7.82	76.2	
N1951-06	DB-14-(20-25)-3-12-22	38	1.19	8.71	9.9	8.28	81.4	
N1951-07	DB-13-(10-15)-3-12-22	39	1.19	8.44	9.63	7.63	76.3	
N1951-08	DB-13-(20-25)-3-12-22	40	1.18	8.70	9.88	7.41	71.6	
N1951-09	DB-8-(15-20)-3-12-22	41	1.16	8.82	9.98	8.52	83.4	
N1951-10	DB-8-(30-35)-3-12-22	42	1.14	8.56	9.7	7.64	75.9	
N1951-11	BLIND-DUPLICATE	43	1.14	8.79	9.93	8.2	80.3	
N1951-11D	BLIND-DUPLICATEDUP	44	1.16	8.79	9.95	8.22	80.3	
N1952-01	WIPE-PIPE	49	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N1952-02	CON-1	50	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N1952-03	CON-2	51	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N1953-01	RO-RBR-200033	45	1.14	8.63	9.77	8.59	86.3	
N1953-02	RO-RBR-200033	46	1.17	8.81	9.98	7.85	75.8	
N1954-01	CL-02-031522	47	1.17	8.54	9.71	8.87	90.2	
N1954-02	CL-02-031522-EPH-2	48	1.19	8.60	9.79	8.92	89.9	
N1956-01	TPK-1	52	1.13	8.62	9.75	8.1	80.9	
N1956-02	TPK-2	53	1.14	8.53	9.67	8.52	86.5	
N1956-03	TPK-3	54	1.15	8.80	9.95	8.35	81.8	
N1956-04	TPK-4	55	1.14	8.57	9.71	8.39	84.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$