



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 3/18/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 03/17/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:25

Out Date: 03/18/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB119197

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1960-15	SB-13-20220316-16.0-17.0	64	1.15	8.64	9.79	6.81	65.5	
N1960-16	SB-13-20220316-18.0-19.0	65	1.18	8.49	9.67	6.1	58.0	
N1968-01	G-2-2-3-4	1	1.14	8.64	9.78	8.72	87.7	
N1968-03	G-2-4-5-6	2	1.18	8.78	9.96	8.87	87.6	
N1968-04	G-2-5-8-9	3	1.14	8.36	9.5	8.18	84.2	
N1968-05	G-2-6-11-12	4	1.14	8.80	9.94	7.42	71.4	
N1968-06	G-2-7-14-15	5	1.13	8.81	9.94	8.5	83.7	
N1968-07	G-2-8-16-17	6	1.18	8.72	9.9	8.08	79.1	
N1968-08	G-2-9-18-19	7	1.19	8.66	9.85	7.93	77.8	
N1968-09	G-2-10-20-21	8	1.19	8.69	9.88	7.3	70.3	
N1968-10	G-2-11-23-24	9	1.12	8.38	9.5	7.07	71.0	
N1968-11	G-2-12-26-27	10	1.16	8.35	9.51	6.94	69.2	
N1968-14	G-3-2-2-3	11	1.17	8.49	9.66	8.96	91.8	
N1968-15	G-3-3-5-6	12	1.19	8.59	9.78	8.61	86.4	
N1970-01	G-3-4-7-8	13	1.17	8.79	9.96	8.95	88.5	
N1970-02	G-3-5-8-9	14	1.19	8.78	9.97	8.96	88.5	
N1970-03	G-3-6-9-10	15	1.18	8.53	9.71	8.56	86.5	
N1970-04	G-3-7-12-13	16	1.15	8.79	9.94	8.26	80.9	
N1970-05	G-3-8-14-15	17	1.17	8.47	9.64	7.68	76.9	
N1970-06	G-3-9-16-17	18	1.18	8.78	9.96	8.55	83.9	
N1970-07	G-3-10-17-18	19	1.13	8.69	9.82	6.67	63.8	
N1970-08	G-3-11-19-20	20	1.19	8.43	9.62	6.58	63.9	
N1970-09	G-3-12-21-22	21	1.17	8.78	9.95	7.82	75.7	
N1970-10	G-3-13-25-26	22	1.19	8.45	9.64	7.88	79.2	
N1970-11	G-4-1-0-1	23	1.18	8.52	9.7	8.54	86.4	
N1970-12	G-4-2-3-4	24	1.19	8.58	9.77	8.52	85.4	
N1970-13	G-4-3-6-7	25	1.19	8.76	9.95	8.88	87.8	



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 3/18/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 03/17/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:25

Out Date: 03/18/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB119197

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1970-15	G-4-4-9-10	26	1.19	8.32	9.51	8.06	82.6	
N1970-16	G-4-5-10-11	27	1.15	8.80	9.95	8.27	80.9	
N1972-01	G-4-6-14-15	28	1.18	8.79	9.97	8.39	82.0	
N1972-02	G-4-7-17-18	29	1.15	8.81	9.96	8.46	83.0	
N1972-03	G-4-8-21-22	30	1.17	8.74	9.91	8.12	79.5	
N1972-05	G-5-2-1-2	31	1.17	8.49	9.66	9.00	92.2	
N1972-06	G-5-3-4-5	32	1.17	8.79	9.96	9.44	94.1	
N1972-07	G-5-4-6-7	33	1.18	8.43	9.61	8.48	86.6	
N1972-08	G-5-5-9-10	34	1.19	8.78	9.97	9.06	89.6	
N1972-09	G-5-6-12-13	35	1.17	8.77	9.94	8.93	88.5	
N1972-10	G-5-7-14-15	36	1.14	8.84	9.98	8.77	86.3	
N1972-11	G-5-8-15-16	37	1.19	8.54	9.73	7.72	76.5	
N1972-12	G-5-9-17-18	38	1.18	8.78	9.96	8.4	82.2	
N1974-01	G-5-11-25-26	39	1.14	8.83	9.97	8.48	83.1	
N1974-02	G-5-12-27-28	40	1.16	8.82	9.98	8.84	87.1	
N1974-03	G-6-1-0-1	41	1.19	8.62	9.81	9.14	92.2	
N1974-05	G-6-3-4-5	42	1.19	8.69	9.88	9.18	91.9	
N1974-06	G-6-4-6-7	43	1.14	8.80	9.94	8.81	87.2	
N1974-08	G-6-5-9-10	44	1.18	8.78	9.96	9.5	94.8	
N1974-09	G-6-6-12-13	45	1.18	8.79	9.97	9.33	92.7	
N1974-10	G-6-7-14-15	46	1.19	8.50	9.69	8.2	82.5	
N1974-12	G-6-9-17-18	47	1.15	8.78	9.93	8.05	78.6	
N1994-01	SB-01-4.5-5.0	48	1.13	8.54	9.67	8.35	84.5	
N1994-02	SB-01-COMP	49	1.18	8.68	9.86	8.52	84.6	
N1994-04	SB-02-4.5-5.0	50	1.18	8.64	9.82	8.46	84.3	
N1994-05	SB-02-COMP	51	1.17	8.33	9.5	7.96	81.5	
N1994-07	SB-03-4.5-5.0	52	1.15	8.83	9.98	8.36	81.7	
N1994-08	SB-03-COMP	53	1.14	8.80	9.94	8.64	85.2	



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 3/18/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 03/17/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:25

Out Date: 03/18/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB119197

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N1994-10	SB-04-4.5-5.0	54	1.16	8.70	9.86	8.53	84.7	
N1994-11	SB-04-COMP	55	1.17	8.50	9.67	8.43	85.4	
N1994-13	SB-05-4.5-5.0	56	1.15	8.82	9.97	9.14	90.6	
N1994-14	SB-05-COMP	57	1.14	8.77	9.91	8.73	86.5	
N1995-01	SB-06-4.5-5.0	58	1.15	8.59	9.74	7.8	77.4	
N1995-02	SB-06-COMP	59	1.15	8.79	9.94	8.5	83.6	
N1995-04	SB-07-4.5-5.0	60	1.18	8.75	9.93	8.54	84.1	
N1995-05	SB-07-COMP	61	1.18	8.34	9.52	8.25	84.8	
N1995-07	SB-08-4.5-5.0	62	1.19	8.53	9.72	8.25	82.8	
N1995-08	SB-06-COMP	63	1.14	8.61	9.75	8.4	84.3	
N2001-01	D-14A	66	1.15	8.60	9.75	8.6	86.6	
N2007-01	SP1	67	1.14	8.84	9.98	8.68	85.3	
N2007-02	SP2	68	1.14	8.52	9.66	8.8	89.9	
N2007-03	SP3	69	1.18	8.80	9.98	8.99	88.8	
N2009-01	13563	70	1.00	1.00	2.00	2.00	100.0	WASTE FUEL
N2011-01	C-A	71	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N2011-02	C-B	72	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N2011-03	C-C	73	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N2011-04	C-D	74	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N2011-05	C-E	75	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$