



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 3/21/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 16:55

In Date: 03/18/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:14

Out Date: 03/19/2022

Weight Check 1.0g: 0.01

Weight Check 10g: 1.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB119240

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2010-01	WC1-220316	1	1.18	8.52	9.7	8.55	86.5	
N2021-01	2279	2	1.00	1.00	2.00	2.00	100.0	oil sample
N2021-03	MHQ-SP	3	1.18	8.57	9.75	8.69	87.6	
N2021-04	MHQ-SP-EPH-2	4	1.18	8.51	9.69	9.06	92.6	
N2025-01	D-154	5	1.00	1.00	2.00	2.00	100.0	GREASE
N2027-01	NB-08-031822	6	1.18	8.60	9.78	9.14	92.6	
N2027-02	NB-08-031822-EPH-2	7	1.12	8.56	9.68	9.09	93.1	
N2030-01	MR-1	8	1.18	8.59	9.77	8.75	88.1	
N2032-01	FDA47D-1-1	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-02	FDA47D-1-2	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-03	FDA47D-2-1	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-04	FDA47D-2-2	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-05	JEB922V-1-1	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-06	JEB922V-1-2	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-07	JEB922V-2-1	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-08	JEB922V-2-2	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-09	KMH239X-1-1	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-10	KMH239X-1-2	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-11	226642202-1-1	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2032-12	226642202-1-2	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2033-01	D-2378	21	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N2033-03	VNJ-237	22	1.15	8.80	9.95	7.76	75.1	
N2033-04	VNJ-237-EPH-2	23	1.17	8.81	9.98	8.31	81.0	
N2039-01	OK-03-031822	24	1.14	8.50	9.64	8.84	90.6	
N2039-02	OK-03-031822-EPH-2	25	1.12	8.44	9.56	8.59	88.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$