



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 04/20/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 04/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119764

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2480-01	WC-16-1A-1-GR	1	1.15	8.40	9.55	8.62	88.9	
N2480-03	WC-16-1A-1	2	1.17	8.54	9.71	8.51	85.9	
N2480-05	WC-16-1A-2-GR	3	1.16	8.54	9.7	8.14	81.7	
N2480-07	WC-16-1A-2	4	1.16	8.81	9.97	8.52	83.5	
N2480-09	WC-16-1A-3-GR	5	1.11	8.50	9.61	7.82	78.9	
N2480-11	WC-16-1A-3	6	1.19	8.49	9.68	8.04	80.7	
N2480-13	WC-16-1A-4-GR	7	1.15	8.70	9.85	8.33	82.5	
N2480-15	WC-16-1A-4	8	1.15	8.66	9.81	8.21	81.5	
N2487-01	UST-1-PB	20	1.15	8.76	9.91	8.04	78.7	
N2487-02	UST-1-N	21	1.15	8.79	9.94	8.62	85.0	
N2487-03	UST-1-E	22	1.17	8.68	9.85	8.05	79.3	
N2487-04	UST-1-S	23	1.19	8.78	9.97	8.5	83.3	
N2487-05	UST-1-W	24	1.17	8.69	9.86	8.36	82.7	
N2489-01	20220369-WEST-CREEK-TUBE-A-1-COMPOSITE	9	1.18	8.54	9.72	4.25	35.9	SLUDGE SAMPLE
N2489-02	20220370-WEST-CREEK-TUBE-A-5-COMPOSITE	10	1.15	8.80	9.95	4.32	36.0	SLUDGE SAMPLE
N2489-03	20220371-WEST-CREEK-TUBE-B-8-COMPOSITE	11	1.15	8.81	9.96	4.05	32.9	SLUDGE SAMPLE
N2489-04	20220372-WEST-CREEK-TUBE-B-14-COMPOSITE	12	1.16	8.74	9.9	4.12	33.9	SLUDGE SAMPLE
N2489-05	20220373-WEST-CREEK-TUBE-B-22-COMPOSITE	13	1.15	8.83	9.98	4.18	34.3	SLUDGE SAMPLE
N2489-06	20220374-WEST-CREEK-TUBE-A-22-COMPOSITE	14	1.14	8.79	9.93	4.69	40.4	SLUDGE SAMPLE
N2491-01	RBR-200046-333	15	1.15	8.83	9.98	8.82	86.9	
N2493-01	72-11944	16	1.15	8.80	9.95	9.05	89.8	
N2493-02	72-11944-EPH-2	17	1.15	8.80	9.95	9.05	89.8	
N2494-01	MH-1S	25	1.14	8.70	9.84	8.87	88.9	
N2494-03	MH-1S-EPH	26	1.14	8.70	9.84	8.87	88.9	
N2494-04	MH-1D	27	1.16	8.59	9.75	9.28	94.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 04/20/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 04/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119764

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2494-06	MH-1D-EPH	28	1.16	8.59	9.75	9.28	94.5	
N2507-01	PL-01-042022	18	1.11	8.68	9.79	9.1	92.1	
N2507-02	PL-01-042022-EPH-2	19	1.11	8.68	9.79	9.1	92.1	
N2508-01	40576-SOL	29	1.15	8.82	9.97	5.59	50.3	
N2510-01	AUD-1596	30	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2510-02	AUD-1597	31	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2510-03	AUD-1598	32	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2510-04	AUD-1599	33	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2510-05	AUD-1600	34	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2510-06	AUD-1602	35	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N2511-01	1526	36	1.14	8.50	9.64	8.21	83.2	
N2511-02	1527	37	1.15	8.40	9.55	9.35	97.6	
N2511-03	1528	38	1.19	8.56	9.75	7.43	72.9	
N2511-04	1529	39	1.16	8.81	9.97	9.71	97.0	
N2511-05	1526-29	40	1.13	8.67	9.8	8.61	86.3	
N2511-05D	1526-29DUP	41	1.15	8.67	9.82	8.63	86.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$