



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 5/13/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 16:55

In Date: 04/27/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:28

Out Date: 04/28/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB119858

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2514-09	WC-16-1A-6	1	1.15	8.53	9.68	8.00	80.3	
N2564-01	PSC-164718	2	1.14	8.52	9.66	9.11	93.5	
N2564-03	PSC-124019	3	1.16	8.51	9.67	7.67	76.5	GEL MATRIX
N2569-01	OR-02-042622	4	1.19	8.72	9.91	7.41	71.3	GEL MATRIX
N2569-02	OR-02-042622-EPH-2	5	1.19	8.72	9.91	7.41	71.3	GEL MATRIX
N2575-12	30899	20	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE,
N2577-01	GT-08-(7.5-10)	6	1.15	8.77	9.92	5.94	54.6	
N2577-02	GT-08-(25-30)	7	1.15	8.81	9.96	8.63	84.9	
N2577-03	GT-06-(8.5-10)	8	1.18	8.70	9.88	8.71	86.6	
N2577-04	GT-06-(25-30)	9	1.14	8.40	9.54	8.48	87.4	
N2577-05	BLIND-DUPLICATE-4-25-22	10	1.15	8.77	9.92	8.16	79.9	
N2577-06	GT-04-(10-15)	11	1.16	8.80	9.96	8.68	85.5	
N2577-07	GT-04-(25-30)	12	1.18	8.75	9.93	9.12	90.7	
N2577-08	GT-11-(15-20)	13	1.13	8.70	9.83	8.32	82.6	
N2577-09	GT-11-(25-30)	14	1.17	8.66	9.83	8.6	85.8	
N2577-10	N2577-07MS	15	1.18	8.75	9.93	9.12	90.7	
N2577-11	N2577-07MSD	16	1.18	8.75	9.93	9.12	90.7	
N2578-01	K22084M-END-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2578-02	K22084M-END-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2588-01	MH-1	17	1.19	8.58	9.77	8.83	89.0	
N2588-03	MH-2	18	1.14	8.61	9.75	8.89	90.0	
N2588-03D	MH-2DUP	19	1.15	8.62	9.77	8.88	89.7	
N2597-01	KMH332L-END-1-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-02	KMH332L-END-1-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-03	KMH332L-END-2-1	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-04	KMH332L-END-2-2	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-05	KMH393V-END-1-1	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-06	KMH393V-END-1-2	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: mohammad
Analyst: jignesh
Date: 5/13/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 16:55
In Date: 04/27/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 04/28/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119858

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2597-07	KMH249C-END-1-1	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2597-08	KMH249C-END-1-2	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2598-01	OS-DRUM-COMP	31	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N2598-02	OS-DRUM-COMP-EPH-2	32	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N2599-01	M-1	33	0.92	9.06	9.98	8.88	87.9	
N2600-01	MH-15S	34	0.91	8.74	9.65	9.24	95.3	
N2600-03	MH-15S-EPH	35	0.9	8.83	9.73	9.33	95.5	
N2600-04	MH-15D	36	0.93	8.78	9.71	9.31	95.4	
N2600-06	MH-15D-EPH	37	0.92	8.99	9.91	9.42	94.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$