



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/5/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 05/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:22
Out Date: 05/05/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2666-01	WASTE	1	0.92	9.05	9.97	7.58	73.6	
N2672-01	01	2	0.93	8.76	9.69	8.07	81.5	
N2672-02	02	3	0.93	8.83	9.76	7.45	73.8	
N2672-03	03	4	0.91	8.77	9.68	7.42	74.2	
N2672-04	04	5	0.94	9.02	9.96	8.65	85.5	
N2672-05	05	6	0.92	9.06	9.98	9.2	91.4	
N2672-06	06	7	0.9	8.83	9.73	7.94	79.7	
N2672-07	07	8	0.94	9.03	9.97	8.00	78.2	
N2672-08	08	9	0.94	8.82	9.76	7.56	75.1	
N2672-09	09	10	0.91	9.00	9.91	8.23	81.3	
N2672-10	10	11	0.95	9.03	9.98	8.59	84.6	
N2676-01	GT-05- (11-15)	12	0.95	8.96	9.91	8.9	88.7	
N2676-02	GT-05- (23-25)	13	0.89	9.04	9.93	7.37	71.7	
N2676-03	GT-10- (10-12)	14	0.93	9.04	9.97	8.26	81.1	
N2676-04	GT-10- (23-25)	15	0.9	8.80	9.7	8.04	81.1	
N2676-04D	GT-10- (23-25) DUP	16	0.93	8.79	9.72	8.06	81.1	
N2680-01	B1-0-2	17	0.92	8.99	9.91	8.51	84.4	
N2680-02	B1-5-7	18	0.93	8.63	9.56	8.29	85.3	
N2680-03	B2-0-2	19	0.95	9.02	9.97	8.63	85.1	
N2680-04	B2-5-7	20	0.94	9.03	9.97	8.55	84.3	
N2680-05	B3-0-2	21	0.91	9.00	9.91	8.25	81.6	
N2680-06	B3-6-8	22	0.96	9.00	9.96	8.66	85.6	
N2680-07	B4-0-2	23	0.94	9.00	9.94	8.71	86.3	
N2680-08	B4-5-7	24	0.96	8.95	9.91	8.27	81.7	
N2680-09	B5-0-2	25	0.88	8.73	9.61	8.95	92.4	
N2680-10	B5-5-7	26	0.92	9.02	9.94	9.16	91.4	
N2684-01	BT-OIL	33	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N2684-03	BT-OILY-DEB	34	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/5/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 05/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:22
Out Date: 05/05/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2684-06	BT-OIL	35	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N2684-07	BT-OILY-DEB	36	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N2685-01	WC-8-19-1-1	52	0.92	8.85	9.77	7.74	77.1	
N2685-02	WC-8-19-1-2	53	0.94	9.00	9.94	8.73	86.6	
N2685-03	WC-8-19-1-3	54	0.9	9.08	9.98	8.24	80.8	
N2685-04	WC-14B-5	55	0.92	8.78	9.7	8.23	83.3	
N2685-06	WC-14B-5	56	0.94	8.81	9.75	8.26	83.1	
N2689-01	OK-01-050422	37	0.92	8.75	9.67	8.5	86.6	
N2689-02	OK-01-050422-EPH-2	38	0.92	8.67	9.59	8.1	82.8	
N2691-01	DEEP-WC-1	27	0.92	8.91	9.83	8.32	83.1	
N2691-02	DEEP-WC-1-EPH-2	28	0.92	8.91	9.83	8.32	83.1	
N2691-03	SHALLOW-WC-1	29	0.94	9.00	9.94	7.77	75.9	
N2691-04	SHALLOW-WC-1-EPH-2	30	0.94	9.00	9.94	7.77	75.9	
N2691-05	SHALLOW-WC-2	31	0.9	8.99	9.89	9.08	91.0	
N2691-06	SHALLOW-WC-2-EPH-2	32	0.9	8.99	9.89	9.08	91.0	
N2692-01	72-12002	39	0.93	8.77	9.7	8.93	91.2	
N2692-02	72-12002-EPH-2	40	0.93	8.77	9.7	8.93	91.2	
N2693-01	15-A	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-02	15-B	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-03	16-A	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-04	16-B	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-05	17-A	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-06	17-B	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-07	18-A	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2693-08	18-B	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
N2694-01	1952-1943-COMP1	49	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N2694-02	1944-1950-COMP2	50	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N2694-03	1949	51	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/5/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 05/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:22
Out Date: 05/05/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB119970

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2703-01	TP-28	57	0.93	8.84	9.77	9.17	93.2	
N2703-03	TP-28-EPH	58	0.89	8.82	9.71	9.17	93.9	
N2703-04	MH-29S	59	0.94	9.04	9.98	9.00	89.2	
N2703-06	MH-29S-EPH	60	0.9	8.97	9.87	8.61	86.0	
N2703-07	MH-29D	61	0.94	8.60	9.54	8.03	82.4	
N2703-09	MH-29D-EPH	62	0.92	9.05	9.97	8.00	78.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$