



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 5/12/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 05/11/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:25

Out Date: 05/12/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120104

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2770-10	M-15848	11	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N2775-01	MH-35S	1	0.92	8.79	9.71	8.99	91.8	
N2775-03	MH-35S-EPH	2	0.93	8.59	9.52	8.67	90.1	
N2775-04	MH-35D	3	0.96	8.66	9.62	8.91	91.8	
N2775-06	MH-35D-EPH	4	0.91	8.60	9.51	8.69	90.5	
N2776-01	CF-280MER-8-15	5	0.93	8.81	9.74	8.72	88.4	
N2776-02	CF-280MER-9-5	6	0.94	8.99	9.93	8.65	85.8	
N2776-03	CF-280MER-10-15	7	0.94	9.01	9.95	9.16	91.2	
N2776-03D	CF-280MER-10-15DUP	8	0.93	9.04	9.97	9.19	91.4	
N2780-01	ELI-2022-19-A	12	1.00	1.00	2.00	2.00	100.0	pilc
N2780-02	ELI-2022-19-B	13	1.00	1.00	2.00	2.00	100.0	pilc
N2780-03	ELI-2022-20-A	14	1.00	1.00	2.00	2.00	100.0	pilc
N2780-04	ELI-2022-21-A	15	1.00	1.00	2.00	2.00	100.0	pilc
N2780-05	ELI-2022-21-B	16	1.00	1.00	2.00	2.00	100.0	pilc
N2780-06	ELI-2022-22-A	17	1.00	1.00	2.00	2.00	100.0	pilc
N2780-07	ELI-2022-22-B	18	1.00	1.00	2.00	2.00	100.0	pilc
N2780-08	ELI-2022-23-A	19	1.00	1.00	2.00	2.00	100.0	pilc
N2780-09	ELI-2022-23-B	20	1.00	1.00	2.00	2.00	100.0	pilc
N2780-10	EFB669Q-1-1	21	1.00	1.00	2.00	2.00	100.0	pilc
N2780-11	EFB669Q-1-2	22	1.00	1.00	2.00	2.00	100.0	pilc
N2780-12	EFB669Q-2-1	23	1.00	1.00	2.00	2.00	100.0	pilc
N2780-13	EFB669Q-2-2	24	1.00	1.00	2.00	2.00	100.0	pilc
N2780-14	EFB6460-1-1	25	1.00	1.00	2.00	2.00	100.0	pilc
N2780-15	EFB6460-1-2	26	1.00	1.00	2.00	2.00	100.0	pilc
N2780-16	KMA964P-1-1	27	1.00	1.00	2.00	2.00	100.0	pilc
N2780-17	KMA964P-1-2	28	1.00	1.00	2.00	2.00	100.0	pilc
N2780-18	BC274414-1-1	29	1.00	1.00	2.00	2.00	100.0	pilc
N2780-19	BC274414-1-2	30	1.00	1.00	2.00	2.00	100.0	pilc



PERCENT SOLID

Supervisor: mohammad
Analyst: jignesh
Date: 5/12/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 05/11/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 05/12/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB120104

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
N2784-01	72-11939	31	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N2792-01	RO-371	9	0.93	8.73	9.66	2.91	22.7	gel matrix
N2792-02	RO-371-EPH-2	10	0.93	8.73	9.66	2.91	22.7	GEL MATRIX
N2795-01	VNJ-239	32	0.9	8.77	9.67	8.76	89.6	
N2795-02	VNJ-239	33	0.9	8.88	9.78	8.79	88.9	
N2802-01	SHALLOW-WC-1	34	0.89	8.77	9.66	8.99	92.4	
N2802-02	SHALLOW-WC-1-EPH-2	35	0.89	8.77	9.66	8.99	92.4	
N2802-03	DEEP-WC-1	36	0.9	8.55	9.45	8.77	92.0	
N2802-04	DEEP-WC-1-EPH-2	37	0.9	8.55	9.45	8.77	92.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$