

**PERCENT SOLID**

Supervisor: mohammad

Analyst: jignesh

Date: 5/27/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 05/26/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 05/27/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120370

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3042-01	JEH62F-1-2	1	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3042-02	JEH62F-2-2	2	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3042-03	KMH493Z-1-2	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3042-04	KMH493Z-2-2	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3042-05	303	5	0.91	9.04	9.95	5.98	56.1	
N3042-06	303-EPH-2	6	0.91	9.04	9.95	5.98	56.1	
N3042-07	302	7	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N3052-01	HEE613R-1-1	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-02	HEE613R-1-2	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-03	HEE613R-2-1	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-04	HEE613R-2-2	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-05	HEE613R-3-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-06	HEE613R-3-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-07	KMA563T-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-08	KMA563T-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-09	61-A	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3052-10	61-B	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3055-01	EDA7916-1-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3055-02	EDA7916-1-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3055-03	ETGI-329	20	0.88	8.86	9.74	8.19	82.5	
N3055-04	ETGI-329-EPH-2	21	0.88	8.86	9.74	8.19	82.5	
N3057-01	72-1193	22	0.91	9.04	9.95	8.87	88.1	
N3057-02	72-1193-EPH-2	23	0.91	9.04	9.95	8.87	88.1	
N3061-01	WASTE	24	0.91	9.06	9.97	8.8	87.1	
N3061-02	VOC	25	0.93	9.03	9.96	8.93	88.6	
N3061-03	1	26	0.93	8.66	9.59	8.63	88.9	
N3061-04	2	27	0.93	9.02	9.95	8.94	88.8	
N3061-05	3	28	0.88	8.91	9.79	7.76	77.2	



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 5/27/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 05/26/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 05/27/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120370

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3061-06	4	29	0.93	8.91	9.84	8.69	87.1	
N3061-07	5	30	0.91	8.66	9.57	8.48	87.4	
N3061-07D	5DUP	31	0.93	8.66	9.59	8.5	87.4	
N3062-01	FS-SOIL-0526	51	0.91	8.94	9.85	9.2	92.7	
N3062-02	FS-SOIL-0526-EPH-2	52	0.93	8.74	9.67	8.96	91.9	
N3067-01	33686	53	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
N3069-01	PR-2-D	34	0.91	9.04	9.95	9.42	94.1	
N3069-02	DISP-2-D	35	0.91	8.82	9.73	9.06	92.4	
N3070-01	SC1	32	0.93	8.84	9.77	8.49	85.5	
N3070-03	IG1	33	0.91	8.92	9.83	8.65	86.8	
N3072-01	RT-4206	36	0.88	9.08	9.96	8.49	83.8	GEL MATRIX
N3072-03	PSC-124035	37	0.93	9.03	9.96	6.52	61.9	GEL MATRIX
N3072-05	PSC-124030	38	0.86	8.99	9.85	7.17	70.2	GEL MATRIX
N3074-01	WASTE	39	0.91	8.87	9.78	8.8	89.0	
N3074-02	VOC	40	0.96	9.02	9.98	8.91	88.1	
N3074-03	1	41	0.91	8.97	9.88	8.89	89.0	
N3074-04	2	42	0.91	9.03	9.94	8.8	87.4	
N3074-05	3	43	0.95	8.99	9.94	9.12	90.9	
N3074-06	4	44	0.89	9.05	9.94	8.89	88.4	
N3074-07	5	45	0.93	8.74	9.67	8.72	89.1	
N3074-08	6	46	0.92	8.99	9.91	8.13	80.2	
N3076-01	RBR-671	54	0.91	9.01	9.92	8.94	89.1	
N3076-02	RBR-671	55	0.89	9.08	9.97	9.12	90.6	
N3077-01	CL-02-052622	56	0.91	8.83	9.74	9.3	95.0	
N3077-02	CL-02-052622-EPH-2	57	0.93	9.01	9.94	9.46	94.7	
N3080-01	WC-4	47	0.91	8.74	9.65	8.74	89.6	
N3080-03	WC-4-EPH-2	48	0.92	8.88	9.8	8.56	86.0	
N3080-04	WC-5	49	0.91	8.93	9.84	8.82	88.6	



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 5/27/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:15

In Date: 05/26/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 05/27/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120370

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3080-06	WC-5-EPH-2	50	0.89	8.63	9.52	8.48	87.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$