



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 6/10/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/09/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:28

Out Date: 06/10/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120590

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3268-01	B-215(36-38)	1	0.91	8.82	9.73	8.85	90.0	
N3269-01	14A-TP-7(7.5-8.0)	2	0.92	8.89	9.81	8.03	80.0	
N3269-02	14A-TP-7-(14.5-15.0)	3	0.93	8.83	9.76	4.18	36.8	
N3269-03	14A-TP-8-(7.5-8.0)	4	0.91	9.04	9.95	8.92	88.6	
N3269-04	14A-TP-8(14.5-15.0)	5	0.91	9.01	9.92	8.32	82.2	
N3270-01	TP-1A	6	0.91	8.59	9.5	8.4	87.2	
N3270-02	TP-1B	7	0.93	8.73	9.66	8.42	85.8	
N3270-03	TP-1C	8	0.91	8.62	9.53	8.49	87.9	
N3270-04	TP-1D	9	0.93	8.71	9.64	8.63	88.4	
N3270-05	TP-1E	10	0.93	8.68	9.61	8.52	87.4	
N3270-06	TP-9A	11	0.92	9.03	9.95	8.88	88.2	
N3270-07	TP-9B	12	0.91	8.73	9.64	8.64	88.5	
N3270-08	TP-9C	13	0.91	9.05	9.96	9.02	89.6	
N3270-09	TP-9D	14	0.93	9.04	9.97	9.15	90.9	
N3270-10	TP-9E	15	0.91	8.92	9.83	8.77	88.1	
N3270-11	TP-11A	16	0.91	8.79	9.7	8.52	86.6	
N3270-12	TP-11B	17	0.95	9.02	9.97	8.81	87.1	
N3270-13	TP-11C	18	0.94	8.70	9.64	8.42	86.0	
N3270-14	TP-11D	19	0.91	8.80	9.71	8.51	86.4	
N3270-15	TP-11E	20	0.93	8.68	9.61	8.46	86.8	
N3270-16	TP-2A	21	0.91	8.92	9.83	8.7	87.3	
N3270-17	TP-3A	22	0.9	8.62	9.52	8.66	90.0	
N3270-18	TP-5A	23	0.91	9.06	9.97	9.11	90.5	
N3271-01	A-1	24	0.92	8.85	9.77	9.11	92.5	
N3271-02	A-2	25	0.92	8.68	9.6	8.55	87.9	
N3271-03	A-3	26	0.93	8.84	9.77	8.9	90.2	
N3271-04	A-4	27	0.91	8.59	9.5	8.69	90.6	
N3271-05	A-5	28	0.91	8.92	9.83	8.73	87.7	



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 6/10/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/09/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:28

Out Date: 06/10/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120590

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3271-06	A-6	29	0.91	8.79	9.7	8.71	88.7	
N3271-06D	A-6DUP	30	0.93	8.79	9.72	8.71	88.5	
N3295-01	OR-03-060922	31	0.96	8.93	9.89	8.85	88.4	
N3295-02	OR-03-060922-EPH-2	32	0.96	8.93	9.89	8.85	88.4	
N3296-01	OK-03-060922	33	0.95	9.02	9.97	9.00	89.2	
N3296-02	OK-03-060922-EPH-2	34	0.95	9.02	9.97	9.00	89.2	
N3299-01	COM-4	37	1.00	1.00	2.00	2.00	100.0	concrete
N3299-07	DEC-7N	38	1.00	1.00	2.00	2.00	100.0	concrete
N3299-08	DEC-7E	39	1.00	1.00	2.00	2.00	100.0	concrete
N3299-09	SC-2-DEC-W	40	1.00	1.00	2.00	2.00	100.0	concrete
N3299-10	SC-2-DEC-N	41	1.00	1.00	2.00	2.00	100.0	concrete
N3299-11	SC-DEC-S	42	1.00	1.00	2.00	2.00	100.0	concrete
N3302-01	I-1	35	0.86	9.06	9.92	8.81	87.7	
N3302-02	I-2	36	0.89	8.66	9.55	8.66	89.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$