

**PERCENT SOLID****Supervisor:** mohammad**Analyst:** jignesh**Date:** 6/14/2022**OVENTEMP IN Celsius(°C):** 107**Time IN:** 17:00**In Date:** 06/13/2022**Weight Check 1.0g:** 1.00**Weight Check 10g:** 10.00**OvenID:** M OVEN-1**OVENTEMP OUT Celsius(°C):** 104**Time OUT:** 08:20**Out Date:** 06/14/2022**Weight Check 1.0g:** 1.00**Weight Check 10g:** 10.00**BalanceID:** M SC-4**Thermometer ID:** %SOLIDS-OVEN

QC:LB120631

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3299-03	SC-15	1	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N3299-04	SC-16	2	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N3299-05	SC-17	3	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N3299-06	SC-18	4	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N3315-01	112-A	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-02	112-B	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-03	BC216222-1-1	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-04	BC216222-1-2	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-05	BC216222-1-3	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-06	FDA375-1-1	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-07	FDA375-1-2	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-08	FDA375-1-3	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-09	BC107651-1-1	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-10	BC107651-1-2	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-11	IRV-A-1-1	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-12	IRV-A-1-2	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-13	IRV-A-2-1	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-14	IRV-A-2-2	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-15	BC226557-1-1	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-16	BC226557-1-2	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-17	BC226557-2-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-18	IRV-B-1-1	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3315-19	IRV-B-1-2	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-01	BC271327-A	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-02	BC271327-B	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-03	BC226596-A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-04	BC226596-B	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-05	BC226569-A	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 6/14/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/13/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:20

Out Date: 06/14/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120631

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3316-06	BC226569-B	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-07	FCA566M-A	30	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-08	FCA566M-B	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3316-09	FCA575L	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3319-01	TR-06-061022	33	0.93	8.65	9.58	9.00	93.3	
N3319-02	TR-06-061022-EPH-2	34	0.93	8.65	9.58	9.00	93.3	
N3327-01	EW-WMR-COMP-1	35	0.91	8.70	9.61	8.65	89.0	
N3327-02	EW-WMR-EPH-2	36	0.9	8.75	9.65	8.55	87.4	
N3327-03	EW-WMR-COMP-2	37	0.91	8.80	9.71	8.59	87.3	
N3327-04	EW-WMR-EPH-4	38	0.91	9.03	9.94	8.76	86.9	
N3327-05	EW-WMR-COMP-3	39	0.93	8.65	9.58	8.29	85.1	
N3327-06	EW-WMR-EPH-6	40	0.91	9.03	9.94	8.71	86.4	
N3328-01	BU-02-061022	41	0.94	8.72	9.66	8.88	91.1	
N3329-01	PALMYRA-RT-73-COMP-1	42	0.9	8.78	9.68	8.63	88.0	
N3329-03	PALMYRA-RT-73-COMP-2	43	0.93	9.03	9.96	9.4	93.8	
N3329-05	PALMYRA-RT-73-COMP-3	44	0.91	9.07	9.98	9.42	93.8	
N3329-07	PALMYRA-RT-73-EPH-1	45	0.93	8.58	9.51	8.92	93.1	
N3329-08	PALMYRA-RT-73-EPH-2	46	0.91	8.82	9.73	8.96	91.3	
N3329-09	PALMYRA-RT-73-EPH-3	47	0.91	8.66	9.57	8.99	93.3	
N3329-10	PALMYRA-RT-73-EPH-4	48	0.91	9.01	9.92	9.28	92.9	
N3329-11	PALMYRA-RT-73-EPH-5	49	0.91	8.92	9.83	9.09	91.7	
N3329-12	PALMYRA-RT-73-EPH-6	50	0.91	9.06	9.97	9.11	90.5	
N3329-13	PALMYRA-RT-73-EPH-7	51	0.92	9.01	9.93	9.37	93.8	
N3329-14	PALMYRA-RT-73-EPH-8	52	0.9	8.88	9.78	9.11	92.5	
N3329-15	PALMYRA-RT-73-EPH-9	53	0.9	9.08	9.98	9.22	91.6	
N3329-16	PALMYRA-RT-73-EPH-10	54	0.9	9.04	9.94	9.16	91.4	
N3329-17	PALMYRA-RT-73-EPH-11	55	0.91	8.75	9.66	9.13	93.9	
N3329-18	PALMYRA-RT-73-EPH-12	56	0.89	9.07	9.96	9.27	92.4	



PERCENT SOLID

Supervisor: mohammad

Analyst: jignesh

Date: 6/14/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/13/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:20

Out Date: 06/14/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120631

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3333-01	D-15-TP-15 (7.5-8.0)	57	0.91	9.04	9.95	8.12	79.8	
N3333-02	D-15-TP-15 (14.5-15.0)	58	0.89	8.94	9.83	6.92	67.4	
N3334-01	WC-14A-1A-7	59	0.88	8.97	9.85	8.91	89.5	
N3334-03	WC-14A-1A-7	60	0.9	8.90	9.8	8.09	80.8	
N3334-05	WC-14A-1A-7-GR	61	0.89	9.05	9.94	8.6	85.2	
N3334-07	WC-14A-1A-8	62	0.89	8.94	9.83	8.82	88.7	
N3334-09	WC-14A-1A-8	63	0.9	9.08	9.98	9.03	89.5	
N3334-11	WC-14A-1A-8-GR	64	0.91	8.98	9.89	8.79	87.8	
N3334-13	WC-14A-1A-11	65	0.91	9.05	9.96	7.94	77.7	
N3334-15	WC-14A-1A-11	66	0.9	9.06	9.96	7.87	76.9	
N3334-17	WC-14A-1A-11-GR	67	0.9	8.82	9.72	8.11	81.7	
N3334-21	WC-14A-1A-5	68	0.9	8.84	9.74	8.72	88.5	
N3334-23	WC-14A-1A-5-GR	69	0.89	8.86	9.75	8.71	88.3	
N3336-01	HD-02-061322	70	0.91	8.79	9.7	8.55	86.9	
N3336-02	HD-02-061322-EPH-2	71	0.91	7.88	8.79	8.55	97.0	
N3337-01	SU-02-061322	72	0.93	8.80	9.73	8.59	87.0	
N3337-02	SU-02-061322-EPH-2	73	0.93	8.80	9.73	8.59	87.0	
N3340-01	CS-1	74	0.91	9.04	9.95	9.38	93.7	
N3340-02	CS-2	75	0.91	8.78	9.69	9.12	93.5	
N3340-03	CS-3	76	0.93	8.82	9.75	8.98	91.3	
N3340-03D	CS-3DUP	77	0.94	8.83	9.77	9.00	91.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$