

**PERCENT SOLID**

Supervisor: mohammad

Analyst: jignesh

Date: 6/15/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 06/14/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:14

Out Date: 06/15/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120661

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3334-19	WC-14A-1A-5	1	0.9	8.79	9.69	8.63	87.9	
N3339-01	CB-224-(27-29)	2	0.91	8.90	9.81	8.75	88.1	
N3339-02	CB-226-(26-28)	3	0.91	8.86	9.77	9.32	94.9	
N3347-01	CL-02-061422	4	0.92	9.06	9.98	8.99	89.1	
N3347-02	CL-02-061422-EPH-2	5	0.92	9.06	9.98	8.99	89.1	
N3348-01	TR-5-061422	17	0.96	8.89	9.85	8.26	82.1	
N3348-02	TR-5-061422-EPH-2	18	0.96	9.02	9.98	8.32	81.6	
N3349-01	40040	19	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N3349-02	40041	20	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N3350-01	KMA831E-1-1	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-02	KMA831E-1-2	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-03	BC260963-1-1	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-04	BC260963-1-2	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-05	BC245029-1-1	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-06	BC245029-1-2	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-07	BC137632-1-1	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3350-08	BC137632-1-2	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3351-01	40582	21	0.95	9.02	9.97	8.58	84.6	
N3351-02	40585	22	0.92	8.87	9.79	9.18	93.1	
N3351-03	40588	23	0.94	8.78	9.72	9.31	95.3	
N3351-04	40591	24	0.96	8.81	9.77	8.86	89.7	
N3351-05	40582-84	25	0.96	8.96	9.92	9.36	93.7	
N3351-06	40585-87	26	0.95	8.99	9.94	9.38	93.8	
N3351-07	40588-90	27	0.97	8.98	9.95	9.31	92.9	
N3351-08	40591-94	28	0.96	9.02	9.98	9.34	92.9	
N3352-01	CLI-2022-FCA575L-1-2	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3353-01	PL-01-061422	15	0.91	8.88	9.79	9.31	94.6	
N3353-02	PL-01-061422-EPH-2	16	0.91	8.88	9.79	9.31	94.6	

