



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/7/2022

OVENTEMP IN Celsius(°C): 108

Time IN: 17:00

In Date: 07/06/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 07/07/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120935

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3582-01	HZA562K-1-1	1	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-02	HZA562K-1-2	2	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-03	BC260356LOK-1-1	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-04	BC260356LOK-1-2	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-05	BC260356LOK-2-1	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-06	BC260356LOK-2-2	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-07	KZZ0264-1-1	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-08	KZZ0264-1-2	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-09	KZZ029M-1-1	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3582-10	KZZ029M-1-2	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3583-01	312-310-313	11	0.92	8.63	9.55	8.29	85.4	
N3584-01	CHRT-26572	12	0.92	9.03	9.95	5.06	45.8	GEL MATRIX
N3584-02	CHRT-26572-EPH-2	13	0.92	9.03	9.95	5.06	45.8	GEL MATRIX
N3584-03	CHRT-2716S	14	0.93	8.75	9.68	9.15	93.9	
N3584-04	CHRT-2716S-EPH-2	15	0.93	8.75	9.68	9.15	93.9	
N3584-05	RB-21199	16	0.91	8.90	9.81	6.83	66.5	
N3584-06	RB-21199-EPH-2	17	0.91	8.90	9.81	6.83	66.5	
N3588-01	KMH4424-1-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-02	KMH4424-1-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-03	BC136902-1-1	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-04	BC136902-1-2	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-05	JEH4585-1-1	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-06	JEH4585-1-2	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-07	BC226482-1-1	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-08	BC226482-1-2	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-09	BC226482-2-1	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-10	BC226482-2-2	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3588-11	BC226480-3-1	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/7/2022

OVENTEMP IN Celsius(°C): 108

Time IN: 17:00

In Date: 07/06/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 07/07/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120935

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3588-12	BC226480-3-2	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3590-01	161	30	1.00	1.00	2.00	2.00	100.0	ool sample
N3590-02	162	31	1.00	1.00	2.00	2.00	100.0	debris
N3590-05	160	46	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N3592-01	KZZ157T-1-1	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-02	KZZ157T-1-2	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-03	KMA670K-1-1	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-04	KMA670K-1-2	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-05	BC259197-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-06	BC259197-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-07	BC259197-2-1	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-08	BC259197-2-2	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-09	BC137632-1-1	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-10	BC137632-1-2	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-11	BC107651-1-1	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-12	BC107651-1-2	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-13	BC245029-1-1	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3592-14	BC245029-1-2	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3593-01	20220483A-1A	47	0.92	8.87	9.79	4.7	42.6	sediment sample
N3593-02	20220483B-1B	48	0.9	8.94	9.84	4.76	43.2	sediment sample
N3593-03	20220483C-1C	49	0.89	8.91	9.8	4.92	45.2	sediment sample
N3593-04	20220484A-2A	50	0.91	8.98	9.89	5.69	53.2	sediment sample
N3593-05	20220484B-2B	51	0.92	8.95	9.87	8.24	81.8	sediment sample
N3593-06	20220484C-2C	52	0.91	8.83	9.74	8.26	83.2	sediment sample
N3593-07	20220485A-3A	53	0.89	8.80	9.69	6.05	58.6	sediment sample
N3593-08	20220485B-3B	54	0.96	8.73	9.69	8.07	81.4	sediment sample
N3593-09	20220485C-3C	55	0.9	8.68	9.58	8.11	83.1	sediment sample
N3593-10	20220486A-4A	56	0.91	9.05	9.96	5.73	53.3	sediment sample



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/7/2022

OVENTEMP IN Celsius(°C): 108

Time IN: 17:00

In Date: 07/06/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 07/07/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120935

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3593-11	20220486B-4B	57	0.86	8.68	9.54	5.63	55.0	sediment sample
N3593-12	20220486C-4C	58	0.91	8.87	9.78	6.22	59.9	sediment sample
N3593-13	20220492A-10A	59	0.94	9.03	9.97	5.52	50.7	sediment sample
N3593-14	20220492B-10B	60	0.91	9.06	9.97	6.9	66.1	sediment sample
N3593-15	20220492C-10C	61	0.91	9.07	9.98	7.71	75.0	sediment sample
N3595-01	20220487A-5A	62	0.93	8.75	9.68	6.72	66.2	sediment sample
N3595-02	20220487B-5B	63	0.93	9.04	9.97	8.09	79.2	sediment sample
N3595-03	20220487C-5C	64	0.88	9.08	9.96	8.98	89.2	sediment sample
N3595-04	20220488A-6A	65	0.91	9.06	9.97	5.89	55.0	sediment sample
N3595-05	20220488B-6B	66	0.91	9.03	9.94	5.23	47.8	sediment sample
N3595-06	20220488C-6C	67	0.91	8.92	9.83	5.31	49.3	sediment sample
N3595-07	20220489A-7A	68	0.93	8.97	9.9	6.39	60.9	sediment sample
N3595-08	20220489B-7B	69	0.91	9.01	9.92	5.39	49.7	sediment sample
N3595-09	20220489C-7C	70	0.95	9.02	9.97	5.35	48.8	sediment sample
N3595-10	20220490A-8A	71	0.88	9.01	9.89	6.45	61.8	sediment sample
N3595-11	20220490B-8B	72	0.93	9.01	9.94	8.13	79.9	sediment sample
N3595-12	20220490C-8C	73	0.93	8.70	9.63	8.25	84.1	sediment sample
N3595-13	20220491A-9A	74	0.93	8.77	9.7	6.61	64.8	sediment sample
N3595-14	20220491B-9B	75	0.86	8.83	9.69	7.38	73.8	sediment sample
N3595-15	20220491C-9C	76	0.92	9.06	9.98	7.98	77.9	sediment sample
N3595-19	N3595-01MS	77	0.93	8.75	9.68	6.72	66.2	sediment sample
N3595-20	N3595-01MSD	78	0.93	8.75	9.68	6.72	66.2	sediment sample
N3597-01	20220506-COMPOSITE-A	79	0.91	8.97	9.88	4.95	45.0	sediment sample
N3597-02	20220507-COMPOSITE-B	80	0.93	9.03	9.96	4.77	42.5	sediment sample
N3597-03	20220508-COMPOSITE-C	81	0.93	9.04	9.97	5.63	52.0	sediment sample
N3597-04	20220509-COMPOSITE-D	82	0.93	8.91	9.84	4.8	43.4	sediment sample
N3597-05	20220510-COMPOSITE-E	83	0.87	9.02	9.89	4.8	43.6	sediment sample
N3597-06	20220511-FIELD-DUPLICATE-COMP	84	0.91	9.06	9.97	5.00	45.1	sediment sample



PERCENT SOLID

Supervisor: mohammad

Analyst: JIGNESH

Date: 7/7/2022

OVENTEMP IN Celsius(°C): 108

Time IN: 17:00

In Date: 07/06/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 07/07/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB120935

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N3602-01	21251	90	1.00	1.00	2.00	2.00	100.0	debris
N3602-03	31530	91	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N3602-05	21245	92	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N3602-07	21244	93	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N3602-09	21247	94	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N3603-01	KMH502E-1-1	85	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3603-02	KMH502E-1-2	86	1.00	1.00	2.00	2.00	100.0	PILC REEL
N3604-01	1975	87	1.00	1.00	2.00	2.00	100.0	debris
N3604-02	1973-1974	88	1.00	1.00	2.00	2.00	100.0	debris
N3604-03	1965-1976	89	0.9	8.69	9.59	9.05	93.8	
N3605-01	OE-2-070622	95	0.92	8.61	9.53	9.06	94.5	
N3605-02	OE-2-070622-EPH-2	96	0.91	8.75	9.66	9.07	93.3	
N3605-03	OE-3-070622	97	0.93	8.77	9.7	9.28	95.2	
N3605-04	OE-3-070622-EPH-2	98	0.91	8.71	9.62	9.09	93.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$