



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/31/2022

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 08/26/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/27/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB121750

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4324-01	BC136902-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4324-02	BC136902-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4324-03	KMH356-1-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4324-04	KMH356-1-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4379-01	SB-01-1.5-2.5	1	0.93	8.90	9.83	8.62	86.4	
N4379-02	SB-01-12-14	2	0.96	8.84	9.8	9.21	93.3	
N4379-03	SB-02-1-2	3	1.08	8.84	9.92	9.42	94.3	
N4379-04	SB-02-10-11	4	1.07	8.75	9.82	9.38	95.0	
N4379-05	SB-03-9-11	5	1.1	8.80	9.9	9.35	93.7	
N4379-06	N4379-05MS	6	1.11	8.79	9.9	9.35	93.7	
N4379-07	N4379-05MSD	7	1.11	8.79	9.9	9.35	93.7	
N4379-08	SB-04-12-13	8	1.08	8.90	9.98	9.66	96.4	
N4379-09	SB-05-2-4	9	1.06	8.89	9.95	9.25	92.1	
N4379-10	SB-05-10-12	10	1.06	8.92	9.98	9.38	93.3	
N4379-11	DUP-X	11	1.07	8.80	9.87	9.47	95.5	
N4379-12	SB-06-9.5-11	12	1.06	8.84	9.9	9.4	94.3	
N4379-13	SB-07-3-5	13	1.07	8.81	9.88	8.95	89.4	
N4379-14	SB-07-12-13	14	1.05	8.93	9.98	8.67	85.3	
N4379-15	SB-08-9-11	15	1.08	8.76	9.84	9.00	90.4	
N4379-16	SB-09-1-2	16	1.07	8.89	9.96	9.4	93.7	
N4379-17	SB-09-10-12	17	1.06	8.89	9.95	9.46	94.5	
N4379-18	SB-10-1.5-2.5	18	1.05	8.92	9.97	9.16	90.9	
N4379-19	SB-10-9-11	19	1.08	8.77	9.85	9.36	94.4	
N4379-20	COMPOSITE-DRUMS	20	1.07	8.89	9.96	8.66	85.4	
N4397-01	VNJ-216	28	0.91	8.67	9.58	9.00	93.3	
N4397-02	VNJ-216	29	0.91	9.07	9.98	9.29	92.4	
N4398-01	ROCK-FINES-QA/QC-COMP	25	0.91	8.86	9.77	9.17	93.2	
N4401-01	59A	30	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/31/2022

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 08/26/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/27/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB121750

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4401-02	59B	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-03	60A	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-04	60B	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-05	61A	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-06	61B	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-07	62A	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-08	63A	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-09	63B	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-10	64A	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-11	64B	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-12	65A	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-13	65B	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-14	66A	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-15	66B	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-16	67A	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-17	67B	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-18	68A	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-19	68B	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4401-20	69A	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4403-01	149TH-STREET-DRUM-5-PP E	26	1.00	1.00	2.00	2.00	100.0	plastic sample
N4403-02	149TH-STREET-DRUM-6-PO LY	27	1.00	1.00	2.00	2.00	100.0	plastic sample
N4414-01	OK-2-082622	50	0.91	8.73	9.64	9.14	94.3	
N4414-02	OK-2-082622-EPH-2	51	0.91	8.73	9.64	9.14	94.3	
N4415-01	OR-02-082622	52	0.96	8.68	9.64	9.22	95.2	
N4415-02	OR-02-082622-EPH-2	53	0.96	8.68	9.64	9.22	95.2	
N4416-01	B	54	0.91	8.93	9.84	8.74	87.7	
N4417-01	F	55	0.91	8.77	9.68	9.04	92.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$