



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/31/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/30/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/31/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB121793

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4433-01	EP-1	1	0.91	8.67	9.58	8.6	88.7	
N4433-02	EP-2	2	0.93	8.63	9.56	9.22	96.1	
N4433-03	EP-3	3	0.91	8.64	9.55	9.1	94.8	
N4433-04	EP-BLIND-DUPLICATE	4	0.93	8.76	9.69	9.22	94.6	
N4433-05	EP-4	5	0.94	8.70	9.64	8.91	91.6	
N4433-06	N4433-05MS	6	0.94	8.70	9.64	8.91	91.6	
N4433-07	N4433-05MSD	7	0.94	8.70	9.64	8.91	91.6	
N4438-01	KTS2-WC-S-TR-01-02-03	8	0.92	8.96	9.88	8.89	89.0	
N4438-02	KTS2-WC-S-MH-01-0-20	9	0.93	9.02	9.95	7.99	78.3	
N4438-03	KTS2-WC-S-MH-02-0-20	10	0.91	9.06	9.97	7.81	76.2	
N4439-01	KTS2-TR-03-0-05	11	0.91	9.03	9.94	9.12	90.9	
N4439-02	KTS2-TR-01-0-5	12	0.95	8.68	9.63	8.69	89.2	
N4439-03	KTS2-MH-02-0-20	13	0.91	8.83	9.74	7.38	73.3	
N4439-04	KTS2-MH-01-0-20	14	0.88	8.68	9.56	7.22	73.0	
N4439-05	KTS2-TR-02-0-12	15	0.88	8.89	9.77	8.38	84.4	
N4443-01	HD-01-083022	16	0.92	8.84	9.76	9.18	93.4	
N4443-02	HD-01-083022-EPH-2	17	0.92	8.84	9.76	9.18	93.4	
N4447-01	AUD-1656-OIL	21	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N4447-03	AUD-1659	22	0.91	8.63	9.54	7.94	81.5	SLUDGE
N4448-01	40080	23	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-02	40081	24	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-03	40082	25	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-04	40083	26	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-05	40084	27	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-06	40085	28	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-07	40086	29	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-08	40087	30	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-09	40088	31	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/31/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/30/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/31/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB121793

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4448-10	40089	32	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-11	40070	33	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-12	40071	34	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-13	40072	35	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-14	40073	36	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-15	40074	37	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-16	40075	38	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-17	40076	39	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-18	40077	40	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-19	40078	41	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-20	40079	42	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-21	33804	43	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-22	33805	44	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4448-23	083022A	45	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4450-01	EO-03-083022	18	0.91	8.60	9.51	8.93	93.3	
N4450-02	EO-03-083022-EPH-2	19	0.91	8.60	9.51	8.93	93.3	
N4451-01	TP1	20	0.93	9.04	9.97	9.00	89.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$