



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 10/5/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 10/04/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 10/05/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB122322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4955-01	LOD-MDL-SOIL-01-QT4-2022	1	1.00	1.00	2.00	2.00	100.0	
N4955-02	LOQ-SOIL-02-QT4-2022	2	1.00	1.00	2.00	2.00	100.0	
N4955-03	MDL-SOIL-03-QT4-2022	3	1.00	1.00	2.00	2.00	100.0	
N4955-04	LOD-MDL-SOIL-04-QT4-2022	4	1.00	1.00	2.00	2.00	100.0	
N4955-05	LOQ-SOIL-05-QT4-2022	5	1.00	1.00	2.00	2.00	100.0	
N4955-06	MDL-SOIL-06-QT4-2022	6	1.00	1.00	2.00	2.00	100.0	
N4955-16	LLOQ-SOIL-QT4-2022	7	1.00	1.00	2.00	2.00	100.0	
N4957-01	SB-17	8	0.92	9.03	9.95	9.92	99.7	
N4957-02	SB-18	9	0.91	8.82	9.73	9.71	99.8	
N4957-03	SB-19	10	0.91	8.75	9.66	9.63	99.7	
N4957-04	SB-20	11	0.92	9.05	9.97	9.94	99.7	
N4957-05	SB-21	12	0.9	8.94	9.84	9.78	99.3	
N4957-06	SB-22	13	0.91	8.92	9.83	9.78	99.4	
N4957-07	SB-23	14	0.91	8.67	9.58	9.56	99.8	
N4957-08	SB-24	15	0.91	9.06	9.97	9.94	99.7	
N4957-09	SB-25	16	0.92	8.67	9.59	9.53	99.3	
N4957-10	SB-26	17	0.91	8.98	9.89	9.85	99.6	
N4957-11	SB-27	18	0.91	8.70	9.61	9.54	99.2	
N4957-12	SB-28	19	0.92	8.79	9.71	9.66	99.4	
N4957-13	SB-29	20	0.91	8.72	9.63	9.56	99.2	
N4957-14	SB-30	21	0.91	8.72	9.63	9.6	99.7	
N4957-15	SB-31	22	0.91	8.69	9.6	9.53	99.2	
N4958-01	33527	23	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-02	33528	24	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-03	33524	25	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-04	33525	26	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-05	33526	27	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 10/5/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 10/04/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 10/05/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB122322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4958-06	34231	28	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-07	34232	29	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-08	34233	30	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-09	34235	31	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-10	34236	32	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-11	34237	33	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-12	34238	34	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-13	34239	35	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-14	34240	36	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-15	34241	37	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-16	34242	38	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-17	34160	39	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-18	34164	40	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-19	34165	41	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-20	34166	42	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-21	34167	43	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-22	34168	44	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4958-23	34169	45	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4960-01	PAT-01-100322	46	0.91	8.70	9.61	9.01	93.1	
N4960-02	PAT-01-100322-EPH-2	47	0.91	8.70	9.61	9.01	93.1	
N4960-03	PAT-02-100322	48	0.91	8.93	9.84	8.93	89.8	
N4960-04	PAT-02-100322-EPH-2	49	0.91	8.93	9.84	8.93	89.8	
N4960-05	PAT-03-100322	50	0.9	8.68	9.58	9.03	93.7	
N4960-06	PAT-03-100322-EPH-2	51	0.9	8.68	9.58	9.03	93.7	
N4961-01	40118	52	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4961-02	40119	53	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4961-03	40121	54	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4961-04	40122	55	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Iwona

Analyst: JIGNESH

Date: 10/5/2022

OVENTEMP IN Celsius(°C): 107

Time IN: 17:00

In Date: 10/04/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:19

Out Date: 10/05/2022

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB122322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4961-05	40123	56	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4962-01	KMH557H-1-1	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-02	KMH557H-1-2	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-03	KMH557H-2-1	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-04	KMH557H-2-2	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-05	KMH557H-3-1	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-06	KMH557H-3-2	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-07	KMA659C-1-1	63	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-08	KMA659C-1-2	64	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-09	KZA317K-1-1	65	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-10	KZA317K-1-2	66	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-11	HIM0524-1-1	67	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-12	HIM0524-1-2	68	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-13	BC226462-1-1	69	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4962-14	BC226462-1-2	70	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4969-03	METER-OIL	71	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N4970-01	PLA-162	72	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4970-02	PLA-163	73	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4972-01	HIH052Y-1-1	74	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4972-02	HIH052Y-1-2	75	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4972-03	BC226462-1-1	76	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4972-04	BC226462-1-2	77	1.00	1.00	2.00	2.00	100.0	PILC REEL
N4985-01	BBC	78	0.91	8.59	9.5	9.36	98.4	
N4986-01	40124	79	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-02	40125	80	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-03	40163	81	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-04	40164	82	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-05	40165	83	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 10/5/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/04/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 10/05/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB122322

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N4986-06	40166	84	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-07	40167	85	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-08	40168	86	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-09	40169	87	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-10	40170	88	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-11	40171	89	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-12	40172	90	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-13	40173	91	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-14	40174	92	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-15	40175	93	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4986-16	40176	94	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N4987-01	AUD-1659	95	0.93	9.05	9.98	8.97	88.8	
N4988-01	AU-5-100422	96	0.93	9.04	9.97	9.00	89.3	
N4995-01	L-1	97	0.91	9.04	9.95	9.14	91.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$