



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/22/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB122973

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N5731-01	LEFT-SIDE	1	0.96	8.98	9.94	8.93	88.8	
N5731-02	RIGHT-SIDE	2	0.91	9.07	9.98	8.8	87.0	
N5733-01	BC274442LOZ-END-1-1	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-02	BC274442LOZ-END-1-2	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-03	BC216342LOZ-END-1-1	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-04	BC216342LOZ-END-1-2	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-05	KMA97IC-END-1-1	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-06	KMA97IC-END-1-2	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-07	LAW-2022-48A	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-08	LAW-2022-49A	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-09	LAW-2022-49B	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-10	LAW-2022-50A	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-11	LAW-2022-50B	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-12	LAW-2022-51A	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-13	LAW-2022-51B	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-14	LAW-2022-52A	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-15	LAW-2022-52B	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-16	LAW-2022-53A	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-17	LAW-2022-53B	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-18	LAW-2022-54A	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-19	LAW-2022-54B	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-20	LAW-2022-55A	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-21	LAW-2022-55B	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-22	LAW-2022-56A	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-23	LAW-2022-56B	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-24	LAW-2022-57A	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-25	LAW-2022-57B	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5733-26	LAW-2022-58A	28	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/22/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB122973

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N5733-27	LAW-2022-58B	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N5734-01	5683	47	1.00	1.00	2.00	2.00	100.0	debris
N5734-06	5617	30	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N5734-07	5686	31	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
N5735-01	TR-04-111822	32	0.92	9.01	9.93	9.14	91.2	
N5735-02	TR-04-111822-EPH-2	33	0.92	9.01	9.93	9.14	91.2	
N5736-01	TRE-1679	34	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N5740-01	MH-12	35	0.89	8.67	9.56	8.39	86.5	
N5740-03	MH-12-EPH	36	0.89	8.67	9.56	8.39	86.5	
N5744-05	SVOC-GPC-BLANK	37	1.00	1.00	2.00	2.00	100.0	
N5744-06	PEST-GPC-BLANK	38	1.00	1.00	2.00	2.00	100.0	
N5744-07	PEST-GPC-BLANK-SPIKE	39	1.00	1.00	2.00	2.00	100.0	
N5744-08	PCB-GPC-BLANK	40	1.00	1.00	2.00	2.00	100.0	
N5744-09	PCB-GPC-BLANK-SPIKE	41	1.00	1.00	2.00	2.00	100.0	
N5744-10	SVOC-GPC2-BLANK	42	1.00	1.00	2.00	2.00	100.0	
N5744-11	PEST-GPC2-BLANK	43	1.00	1.00	2.00	2.00	100.0	
N5744-12	PEST-GPC2-BLANK-SPIKE	44	1.00	1.00	2.00	2.00	100.0	
N5744-13	PCB-GCP2-BLANK	45	1.00	1.00	2.00	2.00	100.0	
N5744-14	PCB-GCP2-BLANK-SPIKE	46	1.00	1.00	2.00	2.00	100.0	
N5749-01	MH-82	48	0.89	9.01	9.9	8.64	86.0	
N5749-03	MH-82-EPH-1	49	0.96	8.97	9.93	8.92	88.7	
N5749-04	MH-82-EPH-2	50	0.92	9.02	9.94	9.3	92.9	
N5749-05	MH-82-EPH-3	51	0.91	8.73	9.64	8.03	81.6	
N5749-06	MH-82-EPH-4	52	0.91	9.06	9.97	8.99	89.2	
N5749-07	MH-82-EPH-5	53	0.91	9.07	9.98	8.45	83.1	
N5749-08	MH-83	54	0.89	8.80	9.69	8.34	84.7	
N5749-10	MH-83-EPH-1	55	0.91	8.83	9.74	8.32	83.9	
N5749-11	MH-83-EPH-2	56	0.9	8.70	9.6	8.14	83.2	



PERCENT SOLID

Supervisor: Iwona
Analyst: JIGNESH
Date: 11/22/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/21/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 11/22/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB122973

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N5749-12	MH-83-EPH-3	57	0.91	8.85	9.76	8.57	86.6	
N5749-13	MH-83-EPH-4	58	0.92	8.66	9.58	8.27	84.9	
N5749-14	MH-83-EPH-5	59	0.89	8.91	9.8	8.5	85.4	
N5749-15	MH-83.5	60	1.13	8.74	9.87	8.09	79.6	
N5749-17	MH-83.5-EPH-1	61	1.12	8.86	9.98	8.6	84.4	
N5749-18	MH-83.5-EPH-2	62	1.12	8.85	9.97	8.49	83.3	
N5749-19	MH-83.5-EPH-3	63	1.15	8.45	9.6	8.14	82.7	
N5749-20	MH-83.5-EPH-4	64	1.12	8.80	9.92	8.41	82.8	
N5749-21	MH-83.5-EPH-5	65	1.12	8.85	9.97	8.17	79.7	
N5749-22	MH-84	66	1.13	8.83	9.96	8.72	86.0	
N5749-24	MH-84-EPH-1	67	1.13	8.84	9.97	8.28	80.9	
N5749-25	MH-84-EPH-2	68	1.13	8.69	9.82	8.23	81.7	
N5749-26	MH-84-EPH-3	69	1.14	8.84	9.98	7.44	71.3	
N5749-27	MH-84-EPH-4	70	1.13	8.85	9.98	8.22	80.1	
N5749-28	MH-84-EPH-5	71	1.13	8.63	9.76	8.17	81.6	
N5751-01	STOCK-PILE	72	1.12	8.43	9.55	8.2	84.0	
N5751-02	STOCK-PILE-EPH-2	73	1.14	8.51	9.65	8.23	83.3	
N5751-03	VNJ-223	74	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N5753-01	72-11986	75	1.12	8.48	9.6	8.32	84.9	
N5753-02	72-11986	76	1.13	8.76	9.89	8.54	84.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$