



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/13/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 12/12/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 12/13/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123208

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N5992-01	AOC2-B4-V1	1	1.11	8.69	9.8	9.77	99.7	
N5992-02	AOC2-B4-S5	2	1.14	8.43	9.57	9.55	99.8	
N5992-03	AOC2-B4-W5	3	1.13	8.57	9.7	9.69	99.9	
N5993-01	AOC1-B5-V1	4	0.92	8.65	9.57	9.54	99.7	
N5993-02	AOC1-N5	5	0.92	8.62	9.54	9.53	99.9	
N5993-03	AOC1-E5	6	1.11	8.39	9.5	9.49	99.9	
N5993-04	AOC1-S5	7	1.17	8.62	9.79	9.75	99.5	
N5993-05	AOC1-W5	8	0.92	8.94	9.86	9.84	99.8	
N5993-06	AOC2-B1-V1	9	0.93	8.81	9.74	9.73	99.9	
N5993-07	AOC2-B1-N5	10	0.91	8.86	9.77	9.76	99.9	
N5993-08	AOC2-B1-W5	11	1.11	8.47	9.58	9.57	99.9	
N5994-01	AOC6-B1	12	0.92	8.92	9.84	9.09	91.6	
N5994-02	AOC6-B2	13	0.92	8.74	9.66	8.94	91.8	
N5994-03	AOC6-B2-E10A	14	1.13	8.55	9.68	8.59	87.3	
N5994-04	AOC6-B4	15	0.89	8.67	9.56	8.63	89.3	
N5994-05	AOC6-B4-W10A	16	0.89	8.86	9.75	8.12	81.6	
N5994-06	AOC6-B5	17	0.96	8.98	9.94	9.31	93.0	
N5994-07	AOC6-B5-N10A	18	0.93	8.88	9.81	9.34	94.7	
N5994-08	AOC6-B2-E10B	19	0.93	9.03	9.96	7.72	75.2	
N5994-09	AOC6-B2-E20A	20	0.93	8.71	9.64	8.99	92.5	
N5994-10	AOC6-B2-E20B	21	1.14	8.54	9.68	7.89	79.0	
N5994-11	AOC6-B4-W10B	22	0.93	8.82	9.75	8.43	85.0	
N5994-12	AOC6-B4-W20A	23	0.92	8.80	9.72	8.72	88.6	
N5994-13	AOC6-B4-W20B	24	0.93	8.60	9.53	7.32	74.3	
N5994-14	AOC6-B5-N10B	25	1.15	8.63	9.78	8.29	82.7	
N5994-15	AOC6-B5-N20A	26	0.89	8.99	9.88	9.37	94.3	
N5994-16	AOC6-B5-N20B	27	0.92	8.73	9.65	8.03	81.4	
N5995-01	1282022-B1-S1	28	1.16	8.74	9.9	9.31	93.2	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/13/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 12/12/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 12/13/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123208

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N6007-01	MOO-CONCRETE-1212	34	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
N6010-01	PCB-120722-01	29	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6010-02	PCB-120722-02	30	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6010-03	PCB-120722-03	31	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6011-01	PCB-12922-01	32	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6012-01	PCB12622-01	33	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6022-01	EO-01-121222	36	0.91	8.79	9.7	8.57	87.1	
N6022-02	EO-01-121222-EPH-2	37	0.91	8.79	9.7	8.57	87.1	
N6023-01	WASTE	35	1.16	8.61	9.77	8.66	87.1	
N6024-01	TP-Q	38	0.91	9.05	9.96	8.9	88.3	
N6024-03	TP-Q-EPH	39	1.15	8.49	9.64	8.37	85.0	
N6024-04	MH-21	40	0.93	8.63	9.56	8.56	88.4	
N6024-06	MH-21-EPH	41	0.93	9.05	9.98	8.86	87.6	
N6025-01	NB-07-121222	42	1.18	8.66	9.84	8.91	89.3	
N6025-02	NB-07-121222-EPH-2	43	1.18	8.66	9.84	8.91	89.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$