



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/19/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/16/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 12/17/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123304

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N6091-01	BU-02-121522	1	0.92	8.97	9.89	8.36	82.9	
N6091-02	BU-03-121522	2	0.96	8.59	9.55	7.86	80.3	
N6092-01	20221254-LOCATION-2	3	0.91	6.35	7.26	3.53	41.3	SLUDGE SAMPLE
N6092-02	20221259-LOCATION-7	4	0.96	9.02	9.98	4.79	42.5	SLUDGE SAMPLE
N6092-03	20221263-LOCATION-11	5	1.14	6.27	7.41	3.2	32.9	SLUDGE SAMPLE
N6092-04	20221267-LOCATION-15	6	0.93	8.61	9.54	4.17	37.6	SLUDGE SAMPLE
N6096-01	34260	7	1.00	1.00	2.00	2.00	100.0	debris
N6096-02	COMP-1	8	0.93	8.62	9.55	8.85	91.9	
N6096-03	33875	9	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N6096-04	33876	10	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N6096-05	33877	11	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N6096-06	33878	12	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
N6097-01	FVSub-TP1-COMP-01	13	0.92	8.83	9.75	8.51	86.0	
N6097-03	FVSub-TP1-COMP-01-EPH-2	14	0.92	8.83	9.75	8.51	86.0	
N6098-01	FVSub-TP3-COMP-04	15	0.92	8.78	9.7	8.26	83.6	
N6098-03	FVSub-TP3-COMP-04-EPH-2	16	0.92	8.78	9.7	8.26	83.6	
N6098-07	FVSub-TP6-COMP-06	17	0.96	9.01	9.97	8.49	83.6	
N6098-09	FVSub-TP6-COMP-06-EPH-2	18	0.96	9.01	9.97	8.49	83.6	
N6098-11	FVSub-TP2-COMP-02	19	1.18	8.78	9.96	9.33	92.8	
N6098-13	FVSub-TP2-COMP-02-EPH-2	20	1.18	8.78	9.96	9.33	92.8	
N6101-01	BC274857-1-1	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-02	BC274857-1-2	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-03	BC274857-2-1	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-04	BC274857-2-2	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-05	BC260212-1-1	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-06	BC260212-1-2	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-07	KZZ178J-1-1	30	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/19/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/16/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 12/17/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123304

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N6101-08	KZZ178J-1-2	31	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-09	KZZ178J-2-1	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-10	KZZ178J-2-2	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-11	BC10651-1-1	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6101-12	BC10651-1-2	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6103-01	TT-210-IDWSO-20221215	21	0.91	8.62	9.53	9.28	97.1	
N6111-01	BC274770-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6111-02	BC274770-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6111-03	BC274770-2-1	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6111-04	BC274770-2-2	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6112-01	OR-03-121622	22	0.93	9.04	9.97	8.78	86.8	
N6112-02	OR-03-121622-EPH-2	23	0.93	9.04	9.97	8.78	86.8	
N6121-01	PCB-121422-01	40	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
N6122-01	BC226734-1-1	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-02	BC226734-1-2	42	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-03	BC271246-1-1	43	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-04	BC271246-1-2	44	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-05	BC271246-1-2	45	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-06	BC271246-2-2	46	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-07	BC274271-1-1	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-08	BC274271-1-2	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-09	BC274271-2-1	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-10	BC274271-2-2	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-11	BC260963-1-1	51	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-12	BC260963-1-2	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-13	BC260963-2-1	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6122-14	BC260963-2-2	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6123-01	65A	55	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/19/2022

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 12/16/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:33
Out Date: 12/17/2022
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123304

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
N6123-02	65B	56	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6123-03	66A	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6123-04	66B	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6124-01	1A	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6124-02	1B	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
N6125-01	B5(0-6)	61	0.93	9.04	9.97	9.08	90.2	
N6126-01	A5(0-6)	62	1.16	8.80	9.96	8.59	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$