

Analysis Method: 9034

ANALYST: Rubina

Parameter: Reactive Sulfide

SUPERVISOR REVIEW BY: Iwona

Run Number: LB123455

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W2941
IODINE SOLUTION .025N 1L	W2924
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB149931BL		1	5.02	50	2.00	0.00	1.92	1.92	0.08	0.00	0.00	12/28/2022	10:30
2	N6199-02		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.06	4.74	12/28/2022	10:33
3	N6199-02DUP		1	5.04	50	2.00	0.00	1.82	1.82	0.18	0.10	7.94	12/28/2022	10:35
4	N6200-02		1	5.02	50	2.00	0.00	1.88	1.88	0.12	0.04	3.19	12/28/2022	13:38
5	N6211-01		1	5.01	50	2.00	0.00	1.84	1.84	0.16	0.08	6.39	12/28/2022	10:40
6	N6212-02		1	5.03	50	2.00	0.00	1.86	1.86	0.14	0.06	4.77	12/28/2022	10:42
7	N6217-01		1	5.03	50	2.00	0.00	1.82	1.82	0.18	0.10	7.95	12/28/2022	10:45
8	N6219-02		1	5.07	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	12/28/2022	10:48
9	N6219-05		1	5.05	50	2.00	0.00	1.82	1.82	0.18	0.10	7.92	12/28/2022	10:50
10	N6220-02		1	5.05	50	2.00	0.00	1.84	1.84	0.16	0.08	6.34	12/28/2022	10:53
11	N6220-05		1	5.03	50	2.00	0.00	1.86	1.86	0.14	0.06	4.77	12/28/2022	10:55
12	N6221-01		1	5.06	50	2.00	0.00	1.88	1.88	0.12	0.04	3.16	12/28/2022	11:00

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume