



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/11/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 01/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123613

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1093-01	TT-214-IDWSO-20230109-1	1	1.12	8.72	9.84	9.63	97.6	
O1093-02	TT-214-IDWSO-20230109-2	2	1.13	8.46	9.59	9.27	96.2	
O1093-03	TT-214-IDWSO-20230109-3	3	1.16	8.68	9.84	9.62	97.5	
O1094-01	STOCKPILE-A	4	1.16	8.72	9.88	8.87	88.4	
O1094-02	STOCKPILE-A-2	5	1.19	8.59	9.78	8.64	86.7	
O1094-03	STOCKPILE-B	6	1.15	8.75	9.9	9.00	89.7	
O1099-01	33892	12	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1099-02	33893	13	1.00	1.00	2.00	2.00	100.0	WIPE SAMPLE
O1100-01	WASTE-117	7	1.16	8.60	9.76	8.62	86.7	
O1102-01	JEC519T-1-1	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-02	JEC519T-1-2	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-03	BC226480-1-1	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-04	BC226480-1-2	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-05	BC226480-2-1	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-06	BC226480-2-2	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-07	BC226467-1-1	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-08	BC226467-1-2	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-09	BC248186-1-1	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-10	BC248186-1-2	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-11	KZA4310-1-1	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-12	KZA4310-1-2	25	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-13	KZA4310-2-1	26	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-14	KZA4310-2-2	27	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-15	KZA4310-3-1	28	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1102-16	KZA4310-3-2	29	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1105-01	TP-6	30	1.14	8.82	9.96	8.06	78.5	
O1105-02	TP-3	31	1.18	8.79	9.97	7.7	74.2	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/11/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/10/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 01/11/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123613

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1106-01	MH-3	8	1.12	8.72	9.84	8.91	89.3	
O1106-03	MH-3-EPH-2	9	1.19	8.62	9.81	8.8	88.3	
O1106-04	TP-B	10	1.19	8.64	9.83	8.86	88.8	
O1106-06	TP-B-EPH-2	11	1.14	8.83	9.97	9.41	93.7	
O1107-01	GZA276A-1-1	32	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-02	GZA276A-1-2	33	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-03	ECA281F-1-1	34	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-04	ECA281F-1-2	35	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-05	KQA110U-1-1	36	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-06	KQA110U-1-2	37	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-07	BC271242-1-1	38	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-08	BC271242-1-2	39	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-09	BC271242-2-1	40	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1107-10	BC271242-2-2	41	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1108-01	2739	42	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1109-01	OR-03-011023	43	1.12	8.70	9.82	9.08	91.5	
O1109-02	OR-03-011023-EPH-2	44	1.12	8.70	9.82	9.08	91.5	
O1111-01	D2856	45	1.00	1.00	2.00	2.00	100.0	debris
O1111-02	D2872	46	1.00	1.00	2.00	2.00	100.0	oily-debris sample
O1111-03	D2873	47	1.00	1.00	2.00	2.00	100.0	OILY DEBRIS
O1112-01	BC260952-1-1	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-02	BC260952-1-2	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-03	KMA946R-1-1	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-04	KMA946R-1-2	51	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-05	BC274770-1-1	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-06	BC274770-1-2	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-07	KMH270P-1-1	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1112-08	KMH270P-1-2	55	1.00	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$