



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/19/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123724

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1173-01	FDH135Y-END-1-1	1	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1173-02	FDH135Y-END-1-2	2	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-01	IRV-2023-1A	3	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-02	IRV-2023-1B	4	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-03	IRV-2023-2A	5	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-04	IRV-2023-2B	6	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-05	IRV-2023-3A	7	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-06	IRV-2023-3B	8	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-07	IRV-2023-4A	9	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-08	IRV-2023-4B	10	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-09	IRV-2023-5A	11	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-10	IRV-2023-5B	12	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-11	IRV-2023-6A	13	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-12	IRV-2023-6B	14	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-13	IRV-2023-7A	15	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-14	IRV-2023-7B	16	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-15	IRV-2023-8A	17	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1174-16	IRV-2023-8B	18	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-01	KEA081H-END-1-1	19	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-02	KEA081H-END-1-2	20	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-03	BC274175LOK-END-1-1	21	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-04	BC274175LOK-END-1-2	22	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-05	BC274175LOK-END-2-1	23	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1177-06	BC274175LOK-END-2-2	24	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1179-01	VNJ-208	25	1.14	8.84	9.98	8.51	83.4	
O1179-02	VNJ-208-EPH-2	26	1.14	8.84	9.98	8.51	83.4	
O1179-03	RB21127	27	1.14	8.47	9.61	8.61	88.2	
O1179-04	RB21127-EPH-2	28	1.14	8.47	9.61	8.61	88.2	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/19/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123724

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1182-01	001	29	1.13	8.52	9.65	9.6	99.4	
O1183-01	001	30	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1183-02	002	31	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1183-03	003	32	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1183-04	004	33	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1183-05	005	34	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1183-06	006	35	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1184-01	APNB-BOT-1	36	1.14	8.77	9.91	9.24	92.4	
O1184-02	APNB-SOUTH	37	1.12	8.67	9.79	9.02	91.1	
O1184-03	APNB-WEST	38	1.13	8.72	9.85	9.41	95.0	
O1184-04	APNB-NORTH	39	1.15	8.81	9.96	8.99	89.0	
O1184-05	APNB-EAST	40	1.16	8.76	9.92	9.45	94.6	
O1185-01	001	41	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1185-02	002	42	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1185-03	003	43	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1185-04	004	44	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1185-05	005	45	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1185-06	006	46	1.00	1.00	2.00	2.00	100.0	OIL SAMPLE
O1188-01	BC260897-END-1-1	47	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1188-02	BC260897-END-1-2	48	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1188-03	BC274857-END-1-1	49	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1188-04	BC274857-END-1-2	50	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1188-05	BC260544LOZ-END-1-1	51	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1188-06	BC260544LOZ-END-1-2	52	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-01	AERIAL-2870-1	53	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-02	AERIAL-2870-2	54	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-03	AERIAL-2870-3	55	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-04	AERIAL-2870-4	56	1.00	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/19/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123724

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1189-05	AERIAL-2870-5	57	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-06	AERIAL-2870-6	58	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-07	AERIAL-2871-1	59	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-08	AERIAL-2871-2	60	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-09	AERIAL-2871-3	61	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-10	AERIAL-2871-4	62	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-11	AERIAL-2871-5	63	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-12	AERIAL-2871-6	64	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-13	AERIAL-2869-1	65	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-14	AERIAL-2869-2	66	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-15	AERIAL-2869-3	67	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-16	AERIAL-2869-4	68	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-17	AERIAL-2869-5	69	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1189-18	AERIAL-2869-6	70	1.00	1.00	2.00	2.00	100.0	PILC REEL
O1190-01	PE-1	71	1.15	8.80	9.95	9.51	95.0	
O1190-02	PE-1D	72	1.15	8.59	9.74	9.19	93.6	
O1190-03	PE-2	73	1.14	8.83	9.97	9.46	94.2	
O1190-04	PE-3	74	1.11	8.87	9.98	9.07	89.7	
O1190-05	PE-4	75	1.14	8.65	9.79	9.44	96.0	
O1190-06	PE-5	76	1.16	8.43	9.59	8.94	92.3	
O1190-07	PE-6	77	1.14	8.84	9.98	9.55	95.1	
O1190-08	PE-7	78	1.14	8.47	9.61	8.9	91.6	
O1190-09	PE-8	79	1.13	8.86	9.99	9.38	93.1	
O1190-10	PE-9	80	1.16	8.64	9.8	9.17	92.7	
O1190-11	PE-10	81	1.16	8.78	9.94	9.39	93.7	
O1193-01	P-011723RB01-BLDG-39-R M-I-C-77	82	1.00	1.00	2.00	2.00	100.0	CAULKING SAMPLE
O1198-01	PL-02-011823	105	1.12	8.52	9.64	8.83	90.5	
O1198-02	PL-02-011823-EPH-2	106	1.12	8.52	9.64	8.83	90.5	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/19/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123724

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1200-01	C-1	83	1.12	8.75	9.87	8.79	87.7	
O1214-01	TJRC-011723-B3-0-10	107	1.15	8.80	9.95	9.24	91.9	
O1214-02	TJRC-011723-B3-0-10	108	1.12	8.49	9.61	8.96	92.3	
O1214-04	TJRC-011723-B3-10-20	109	1.13	8.56	9.69	7.92	79.3	
O1214-05	TJRC-011723-B3-10-20	110	1.14	8.84	9.98	9.03	89.3	
O1214-07	TJRC-011723-B3-0-10-DU P	111	1.14	8.47	9.61	8.95	92.2	
O1214-08	TJRC-011723-B3-0-10-DU P	112	1.14	8.81	9.95	9.39	93.6	
O1214-10	TJRC-011823-B3-20-30	113	1.14	8.84	9.98	8.69	85.4	
O1214-11	TJRC-011823-B3-20-30	114	1.15	8.82	9.97	8.81	86.8	
O1214-13	TJRC-011823-B3-30-40	115	1.13	8.86	9.99	8.5	83.2	
O1214-14	TJRC-011823-B3-30-40	116	1.13	8.63	9.76	8.16	81.5	
O1215-01	BU-02-011823	117	1.14	8.50	9.64	8.88	91.1	
O1217-01	BORING-4	84	1.12	8.71	9.83	8.46	84.3	
O1217-03	BORING-4-EPH-2	85	1.13	8.64	9.77	8.54	85.8	
O1217-14	BORING-2	86	1.15	8.81	9.96	8.32	81.4	
O1217-16	BORING-2-EPH-2	87	1.13	8.79	9.92	8.72	86.3	
O1218-01	WC-2R (0-5)	88	1.14	8.80	9.94	9.08	90.2	
O1220-01	MH-8	89	1.15	8.81	9.96	9.12	90.5	
O1220-03	MH-8-EPH-2	90	1.13	8.74	9.87	8.92	89.1	
O1220-04	TP-G	91	1.17	8.81	9.98	9.00	88.9	
O1220-06	TP-G-EPH-2	92	1.1	8.56	9.66	8.88	90.9	
O1222-01	T-1	93	1.15	8.46	9.61	9.17	94.8	
O1222-02	T-2	94	1.14	8.38	9.52	9.18	95.9	
O1222-03	T-3	95	1.14	8.51	9.65	9.13	93.9	
O1222-04	T-4	96	1.13	8.39	9.52	9.19	96.1	
O1222-05	T-5	97	1.14	8.61	9.75	9.45	96.5	
O1222-06	T-6	98	1.15	8.38	9.53	9.36	98.0	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/19/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 01/18/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 01/19/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123724

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
O1222-07	T-7	99	1.13	8.85	9.98	9.85	98.5	
O1222-08	T-8	100	1.14	8.41	9.55	9.1	94.6	
O1222-09	T-9	101	1.14	8.71	9.85	9.48	95.8	
O1222-10	T-10	102	1.15	8.45	9.6	9.24	95.7	
O1222-11	PIPE-11	103	1.14	8.62	9.76	9.42	96.1	
O1222-12	FILL-12	104	1.15	8.72	9.87	7.91	77.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$