



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/26/2023

OVENTEMP IN Celsius(°C): 107
Time IN: 15:10
In Date: 01/25/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 01/26/2023
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB123798

Table with 9 columns: Lab ID, Client SampleID, Dish #, Dish Wt(g) (A), Sample Wt(g), Dish + Sample Wt(g) (B), Dish+Dry Sample Wt(g) (C), % Solid, Comments. It contains 22 rows of data for samples O1296-01 through O1296-22.

% Solid = (C-A) \* 100 / (B-A)